

Anesthesia

For administering local anesthetic prior to ophthalmic surgery
Do not recap needle after use. Discard all needles according to recognized procedures.


Atkinson Tip Anesthesia Needles Single facet with a 22° angle point.



I425 Peribulbar Needle - Atkinson
25G x 7/8 in (.50 x 22mm)



I427 Peribulbar Needle - Atkinson
27G x 7/8 in (.40 x 22mm)



I023 Anesthesia Needle - Atkinson
23G x 1 1/4 in (.60 x 32mm)
May be used for peribulbar or retrobulbar administration




I025 Anesthesia Needle - Atkinson
25G x 1 1/4 in (.50 x 32mm)
May be used for peribulbar or retrobulbar administration




I027 Anesthesia Needle - Atkinson
27G x 1 1/4 in (.40 x 32mm)
May be used for peribulbar or retrobulbar administration



I123 Retrobulbar Needle - Atkinson
23G x 1 1/2 in (.60 x 38mm)



I125 Retrobulbar Needle - Atkinson
25G x 1 1/2 in (.50 x 38mm)




I127 Retrobulbar Needle - Atkinson
27G x 1 1/2 in (.40 x 38mm)

Anesthesia

Sharp Tip Needles Trifacet with a 12° angle point.



I225 Retrobulbar Needle - Sharp
25G x 1 1/2 in (.50 x 38mm)




I230 Facial Infiltration Needle - Sharp
30G x 5/8 in (.30 x 16mm)




I525 Scleral Track Needle
25G x 3/8 in (.50 x 16 mm)


Sub-Tenon's Anesthesia For administering anesthetic agents into the posterior sub-Tenon's space.



I320 Sub-Tenon's Cannula
20G x 1 in (.90 x 25mm)



I327 Sub-Tenon's Cannula
27G x 7/8 in (.40 x 22mm)



I330 Sub-Tenon's Cannula
30G x 7/8 in (.30 x 22mm)