



# Anterior Chamber Maintenance Cannulas


## Anterior Chamber Maintainers

Self-retaining cannulas designed for the continual infusion of fluid to maintain the AC during cataract surgery or anterior vitrectomy.

Code	Gauge and size	Description	Tip
			
M2658	23G x 9/64"	0.64 x 3.5mm Anterior Chamber Maintainer. 150mm silicone tubing. Annular Rings allow cannula to be self retaining (Lewicky)	
M2657	20G x 9/64"	0.9 x 3.5mm Anterior Chamber Maintainer. 150mm silicone tubing. Annular Rings allow cannula to be self retaining (Lewicky)	


## Anterior Chamber Maintenance Cannulas

For maintaining and forming the anterior chamber during surgery by injecting or removing viscoelastic materials, BSS or air.

Code	Gauge and size	Description	Tip
			
M2639	25G x 1"	0.5 x 25mm Anterior Chamber Cannula. Angled 35°, 7.5mm from tip (Rycroft)	
M2270	30G x 7/8"	0.3 x 22mm Anterior Chamber Cannula. Angled 35°, 4mm from tip (Rycroft)	
M2270A	30G x 7/8"	0.3 x 22mm Anterior Chamber Cannula. Angled 35°, 8mm from tip (Rycroft)	
M2273	30G x 7/8"	0.3 x 22mm Anterior Chamber Cannula. Angled 35°, 5.5mm from tip (Rycroft)	
M2273A	27G x 7/8"	0.4 x 22mm Anterior Chamber Cannula. Angled 35°, 4.5mm from tip (Rycroft)	
M2279	27G x 7/8"	0.4 x 22mm Anterior Chamber Cannula. Angled 45°, 11mm from tip (Rycroft)	
M2275	25G x 7/8"	0.5 x 22mm Anterior Chamber Cannula. Angled 35°, 5mm from tip (Rycroft)	
M2273G	23G x 7/8"	0.64 x 22mm Anterior Chamber Cannula. Angled 35°, 5.5mm from tip (Rycroft)	
M2642	21G x 7/8"	0.81 x 22mm Anterior Chamber Cannula. Angled 35°, 7.5mm from tip (Rycroft)	
M2273B	20G x 7/8"	0.9 x 22mm Anterior Chamber Cannula. Angled 45°, 7.5mm from tip (Rycroft)	

## Transzonular Injection Cannula

For the intraoperative injection of steroid-antibiotic combinations.

Code	Gauge and size	Description	Tip
			
M2280	27G x 7/8"	0.4 x 22mm Transzonular Injection Cannula. Angled 33°, 2mm from tip.	