



Backflush/ Extrusion Handpieces

Backflush..... 40
Extrusion..... 42



Backflush/Extrusion Handpieces

Backflush

SUJ08

20GA Backflush/Extrusion Handpiece

for active & passive aspiration with blunt tip



20GA handpiece with blunt tip, grey colour coding ring. Delivered with silicone tube with male and female luer connector either for connection to the machine or for connection to the handpiece without use of the machine to facilitate flow outside of immediate surgical field. Male luer slip coupler for connection to machine tubing. Non-vented male luer cap, optional for use as backflush handpiece.



SUJ09

20GA Backflush/Extrusion Handpiece

for active & passive aspiration with soft tip brush



20GA handpiece with soft tip brush, grey colour coding ring. Delivered with silicone tube with male and female luer connector either for connection to the machine or for connection to the handpiece without use of the machine to facilitate flow outside of immediate surgical field. Male luer slip coupler for connection to machine tubing. Non-vented male luer cap, optional for use as backflush handpiece.



SUJ10

20GA Backflush/Extrusion Handpiece

for active & passive aspiration with soft tip



20GA handpiece with soft tip, grey colour coding ring. Delivered with silicone tube with male and female luer connector either for connection to the machine or for connection to the handpiece without use of the machine to facilitate flow outside of immediate surgical field. Male luer slip coupler for connection to machine tubing. Non-vented male luer cap, optional for use as backflush handpiece.



SUJ11

23GA Backflush/Extrusion Handpiece

for active & passive aspiration with blunt tip



23GA handpiece with blunt tip, green colour coding ring. Delivered with silicone tube with male and female luer connector either for connection to the machine or for connection to the handpiece without use of the machine to facilitate flow outside of immediate surgical field. Male luer slip coupler for connection to machine tubing. Non-vented male luer cap, optional for use as backflush handpiece.



Backflush/Extrusion Handpieces

Backflush

SUJ12

23GA Backflush/Extrusion Handpiece

for active & passive aspiration with soft tip brush



23GA handpiece with soft tip brush, green colour coding ring. Delivered with silicone tube with male and female luer connector either for connection to the machine or for connection to the handpiece without use of the machine to facilitate flow outside of immediate surgical field. Male luer slip coupler for connection to machine tubing. Non-vented male luer cap, optional for use as backflush handpiece.



SUJ13

23GA Backflush/Extrusion Handpiece

for active & passive aspiration with soft tip



23GA handpiece with soft tip, green colour coding ring. Delivered with silicone tube with male and female luer connector either for connection to the machine or for connection to the handpiece without use of the machine to facilitate flow outside of immediate surgical field. Male luer slip coupler for connection to machine tubing. Non-vented male luer cap, optional for use as backflush handpiece.



SUJ14

25GA Backflush/Extrusion Handpiece

for active & passive aspiration with blunt tip



25GA handpiece with blunt tip, blue colour coding ring. Delivered with silicone tube with male and female luer connector either for connection to the machine or for connection to the handpiece without use of the machine to facilitate flow outside of immediate surgical field. Male luer slip coupler for connection to machine tubing. Non-vented male luer cap, optional for use as backflush handpiece.



SUJ17

27GA Backflush/Extrusion Handpiece

for active & passive aspiration with blunt tip



27GA handpiece with blunt tip, blue colour coding ring. Delivered with silicone tube with male and female luer connector either for connection to the machine or for connection to the handpiece without use of the machine to facilitate flow outside of immediate surgical field. Male luer slip coupler for connection to machine tubing. Non-vented male luer cap, optional for use as backflush handpiece.



Backflush/Extrusion Handpieces

Extrusion

SUJ02

25GA Passive Aspiration Handpiece



Passive Aspiration Handpiece 25 Gauge ready for use with 1.5 mm soft tip cannula. For passive aspiration of fluid from the eye during vitreoretinal procedures. Soft flexible tip on the cannula designed atraumatic entry through retinal or macular tears or holes. Single-use, packaged sterile, 12/box.



SUJ03

23GA Passive Aspiration Handpiece



Passive Aspiration Handpiece 23 Gauge ready for use with 1.5 mm soft tip cannula. For passive aspiration of fluid from the eye during vitreoretinal procedures. Soft flexible tip on the cannula designed for atraumatic entry through retinal or macular tears or holes. Single-use, packaged sterile, 12/box.



SUJ06

23GA Passive Aspiration Handpiece



Passive Aspiration Handpiece 23 Gauge ready for use with blunt cannula. For passive aspiration of fluid from the eye during vitreoretinal procedures. Soft flexible tip on the cannula designed for atraumatic entry through retinal or macular tears or holes. Single-use, packaged sterile, 12/box.



SUJ07

25GA Passive Aspiration Handpiece



Passive Aspiration Handpiece 25 Gauge ready for use with blunt cannula. For passive aspiration of fluid from the eye during vitreoretinal procedures. Soft flexible tip on the cannula designed for atraumatic entry through retinal or macular tears or holes. Single-use, packaged sterile, 12/box.

