

ACCESSORIES

DESIGNED, ENGINEERED AND MANUFACTURED SINCE 1980 
MADE IN GERMANY

www.brumaba.de


BRUMABA
OPERATING TABLE SYSTEMS



ACCESSORIES


FEEL LIKE **ON CLOUD NINE**



BRUMABA

OPERATING TABLE SYSTEMS

BRUMABA is a second-generation, family-owned company founded in 1980. **BRUMABA** develops, manufactures, and globally distributes operating tables, operating chairs, and patient transporters and appropriate accessories for specialty departments in hospitals, clinics, and doctor's offices.



BRUMABA'S CORE COMPETENCE:
MOBILE, SPECIALIZED OPERATING TABLES.

BRUMABA is a professional and competent address as well as adviser for all matters concerning ergonomics and positioning on operating tables. We offer our customers products for specialized departments such as eye and oral surgery, OMS, and plastic surgery, dermatology, as well as for general operations, both for in-patients and out-patients.

BRUMABA'S CORE BENEFIT:
OPERATING TABLES FOR OPERATING WITHOUT
TIRING AND OPTIMUM PATIENT POSITIONING.

This complete catalog of accessories contains appropriate accessories for expanding and customizing your entire specialty operating table system according to medical discipline.



CONTENTS

POSITIONING | HEAD

Head plate, large, 400 x 300 mm _____	8
Head plate, small, 320 x 240 mm _____	8
Horseshoe-shaped headrest, (included as standard) _____	8
Horseshoe-shaped headrest, small _____	9
Horseshoe-shaped headrest with long clamping piece, 110 mm _____	9
Silicone head restraining belt for fastening on the headrest _____	9
Hair transplant support _____	10
Neck cushion, viscoelastic, black _____	10
Head plate "rubins pillow" with ball-and-socket joint _____	10
Headrest, detachable, GENIUS (surcharge) _____	11
Headrest, detachable _____	11
Radiolucent headrest, detachable (surcharge) _____	11
External mounting with head plate, radiolucent _____	12
Microsurgery ring, with integrated mounting bracket _____	12

POSITIONING | ARM / SHOULDER / TORSO / LEG

Microsurgery wrist rest (GENIUS EYE) _____	13
Lower arm support for ophthalmology _____	13
Tongue holder, extended, with flat plate rail _____	13
Armrest OP with fastening belt _____	14
Arm support lateral with fastening belt _____	14
Patient handle, hinged _____	14
Hand surgery table with mounting bracket special "XL" _____	15
Hand table telescopic support, single _____	15
Hand surgery table with double telescopic support, with mounting bracket special "XL" _____	15
Back section cushion with arm support, visco _____	16
Horseshoe-shaped headrest for shoulder arthroscopy with head restraining belt _____	16
Back section extension for shoulder arthroscopy _____	16
Side guard, hinged _____	17
Lateral positioning support, deep _____	17
Lateral positioning support, dorsal (back) _____	17
Body belt splitted, silicone with 2 mounting gliders _____	18
Body belt, silicone with 2 mounting gliders _____	18

POSITIONING | LEG / KNEE

Leg holder, OP _____	19
Mounting adapter with 90° direction change _____	19
Leg holder with 2 ball-and-socket joints _____	19
Joint support _____	20
Arthroscopic leg support _____	20
Padded insert for arthroscopic surgery leg support _____	20
Standing plate (1 pair) _____	20

POSITIONING PADDING

Lordosis support padding _____	22
Half roll padding, small _____	22
Half roll padding, large _____	22
Stomach positioning padding, trapezoid _____	23
Mattress, viscoelastic _____	23
Mattress, viscoelastic, short, with gynecological cutout _____	23

OPERATING TABLE | EXTENSIONS

Back section extension, narrow (OP) _____	24
Back section extension, narrow (GENIUS) _____	24
Back section extension, wide (OP) _____	24
Back section extension, wide (GENIUS) _____	24
Back section extension, extra wide _____	25
Back section extension (GENIUS EYE / GENIUS EYE PRO) _____	25
Foot section extension, wide (OP) _____	25
Foot section extension, wide (GENIUS) _____	25
Side widener / foot section extension, narrow (GENIUS) _____	26
Side extension / foot section extension, narrow (OP) _____	26
Side extension, narrow _____	26
Side extension, radiolucent _____	26
Lower leg plate, combined (VARIUS) _____	27
Lower leg plate, combined, extended (VARIUS) _____	27
Lower leg plates, extended, 1 pair (VARIUS) _____	27

INSTRUMENT TABLES

Tray for motor, swifeling, on column _____	28
Instrument table, mobile, 600 x 400 mm _____	28
Instrument table, mobile, 700 x 500 mm _____	28
Instrument table on the column, 600 x 400 mm _____	29
Instrument table on the column, 700 x 500 mm _____	29
Tray for instrument table, aluminium _____	29



ANESTHESIA / INFUSION STANDS

Anesthetic frame bendable, foldable _____	30
Anesthetic frame, bendable - with oxygen supply _____	30
Anesthesiology tube holder _____	31
Anesthesiology tube holder with remote control holder _____	31
Anesthetic frame, universal _____	32
Infusion stand with ball-and-socket joint _____	32
Infusion stand on the column, hinged _____	33
Infusion stand - swivel arm, extended (360 mm) _____	33
Infusion stand mounted on instrument table _____	33

CONTROL ELEMENTS

Hand remote control type IV _____	34
Remote control holder on infusion stand _____	34
Remote control holder, flexible _____	35
Remote control holder, short _____	35
Magnetic plate for remote control _____	35
Foot control, type IV _____	36
Foot control, type IV without memory _____	36
Platform for foot switch type IV - 4-EASY DRIVE _____	37
Foot control, type IV fixed on undercarriage _____	37
Foot control, type IV EYE, without memory _____	37

MOUNTING | BRACKETS AND ADAPTERS

Mounting bracket _____	38
Mounting bracket special "XL" _____	38
Mounting glider with clamping _____	39
Mounting adapter, 25 x 10 mm _____	39
Mounting adapter, strengthened, 290 mm _____	39

MOUNTING | MOUNTING RAILS

Mounting rail, 250 mm head section (GENIUS) _____	40
Mounting rail, 340 mm foot section (GENIUS) _____	40
Mounting rail, 250 mm, massiv seat section (GENIUS) _____	40
Mounting rail, 340 mm, solid back and foot section (GENIUS) _____	41
Mounting rail, 125 mm head section (GENIUS EYE, -EYE PRO) _____	41
Mounting rail set (VENUS) _____	41

TRANSPORT HANDLES

Transport handles, foldaway, head end (1 pair) (GENIUS) _____	42
Transport handles foldaway, head end, (1 pair) (FOCUS, PRIMUS, VARIUS) _____	42
Transport handles, hinged on foot side (1 pair) _____	42

HYGIENE / PRESERVATION / CARE

EasyCover (150 pieces - disposable, unsterile) _____	43
Paper towel holder _____	43
Paper towel holder, 550 mm (VENUS) _____	43
Protective cover for foot section, visco _____	44
Protective cover for foot section, visco, (VARIUS) (1 pair) _____	44
Service kit _____	44

MAGNIFYING LAMPS

Magnifying lamp with swivel arm _____	45
Magnifying lamp with swivel arm (VENUS) _____	45


ACCESSORIES ORDER

Wall mounting rail _____	46
Gliders for wall mounting rail / accessories trolley _____	46
Wall mount for hand surgery table _____	46
Accessories column _____	46
Accessories column - storage plate _____	47
Accessories trolley _____	47
Transport trolley for batteries _____	47



HEAD PLATE, LARGE, 400 x 300 mm

ORDER NUMBER: Z1.003



PADDING DIMENSIONS (L x W x H):
400 x 300 x 50 mm

CLAMPING PIECE HEIGHT: **110 mm**



- ▮ For head positioning, lateral positioning
- ▮ Abdominal positioning in conscious patients
- ▮ Ball joint for seamless inclination, height-adjustable

Suitable for:

- » FOCUS » GENIUS
- » PRIMUS » GENIUS EYE
- » VARIUS » GENIUS EYE PRO

HEAD PLATE, SMALL, 320 x 240 mm

ORDER NUMBER: Z1.109



PADDING DIMENSIONS (L x W x H):
320 x 240 x 50 mm

CLAMPING PIECE HEIGHT: **110 mm**


- ▮ For head positioning, lateral positioning
- ▮ Abdominal positioning in conscious patients
- ▮ Ball joint for seamless inclination, height-adjustable

Suitable for:

- » FOCUS » GENIUS
- » PRIMUS » GENIUS EYE
- » VARIUS » GENIUS EYE PRO

HORSESHOE-SHAPED HEADREST (INCLUDED AS STANDARD)

ORDER NUMBER: Z1.043



PADDING DIMENSIONS (L x W x H):
250 x 240 x 110 mm

CLAMPING PIECE HEIGHT: **65 mm**

- ▮ For a defined and stable head position
- ▮ Ball-and-socket joint, height adjustable

Suitable for:

- » FOCUS » GENIUS
- » PRIMUS » GENIUS EYE
- » VARIUS » GENIUS EYE PRO

HORSESHOE-SHAPED HEADREST, SMALL

ORDER NUMBER: Z1.043A



PADDING DIMENSIONS (L x W x H):
200 x 210 x 80 mm

CLAMPING PIECE HEIGHT: **65 mm**

- ▮ For good access to the ears
- ▮ For head positioning in children
- ▮ Ball-and-socket hinge, height adjustable,

Suitable for:

- » FOCUS » GENIUS
- » PRIMUS » GENIUS EYE
- » VARIUS » GENIUS EYE PRO

HORSESHOE-SHAPED HEADREST WITH LONG CLAMPING PIECE, 110 mm

ORDER NUMBER: Z1.249



CLAMPING PIECE HEIGHT: **110 mm**

- ▮ For Bekhterev's disease patients
- ▮ Defined and stable head position
- ▮ Ball-and-socket hinge, height adjustable

Suitable for:

- » FOCUS » GENIUS
- » PRIMUS » GENIUS EYE
- » VARIUS » GENIUS EYE PRO

SILICONE HEAD RESTRAINING BELT FOR FASTENING ON THE HEADREST

ORDER NUMBER: Z1.189



BELT SIZE (L x W):
900 x 50 mm

- ▮ Hygienic silicone belt for the fixation of the patient's head

Suitable for:

- » FOCUS » GENIUS
- » PRIMUS » GENIUS EYE
- » VARIUS » GENIUS EYE PRO

HAIR TRANSPLANT SUPPORT

ORDER NUMBER: Z1.032



PADDING DIMENSIONS (Ø x L):
90 x 230 mm

- ▮ For optimum access to the upper and rear head area
- ▮ Ball-and-socket hinge, height adjustable

Suitable for:

- » FOCUS » GENIUS
- » PRIMUS » GENIUS EYE
- » VARIUS » GENIUS EYE PRO

NECK CUSHION, VISCOELASTIC, BLACK

ORDER NUMBER: Z1.250



PADDING DIMENSIONS (Ø x L):
200 x 200 mm

CLAMPING PIECE HEIGHT: **65 mm**

- ▮ For extra soft neck positioning
- ▮ Ball-and-socket hinge, height adjustable

Suitable for:

- » FOCUS » GENIUS
- » PRIMUS » GENIUS EYE
- » VARIUS » GENIUS EYE PRO

HEAD PLATE “RUBINS PILLOW” WITH BALL-AND-SOCKET JOINT

ORDER NUMBER: X.000184



PADDING DIMENSIONS (L x W x H):
470 x 250 x 100 mm

CLAMPING PIECE HEIGHT: **110 mm**

- ▮ For defined head positioning
- ▮ Armrest for the surgeon during eye surgery in 12 o'clock position
- ▮ Ball-and-socket hinge, height adjustable,

Suitable for:

- » FOCUS » GENIUS
- » PRIMUS » GENIUS EYE
- » VARIUS » GENIUS EYE PRO

HEADREST, DETACHABLE GENIUS (SURCHARGE)

ORDER NUMBER: X.000257



- ! For a quick change of the head section
- ! Good access for the anesthetist, e.g. during knee or gynecological surgery

Suitable for:

- » GENIUS

HEADREST OP, DETACHABLE

ORDER NUMBER: X.000209



- ! For a quick change of the head section
- ! Good access for the anesthetist, e.g. during knee or gynecological surgery
- ! Solid stainless steel construction

Suitable for:

- » FOCUS
- » PRIMUS
- » VARIUS

RADIOLUCENT HEADREST, DETACHABLE (SURCHARGE)

ORDER NUMBER: X.000235




OPTIONAL: PADDING FOR HEAD SECTION RADIOLUCENT (AS SHOWN)

ORDER NUMBER: E.002553

- ! For a quick change of the head section
- ! Good for intra-operative X-ray, e.g. ENT, OMS

Suitable for:

- » FOCUS
- » PRIMUS
- » GENIUS
- » VARIUS



EXTERNAL MOUNTING WITH HEAD PLATE, RADIOLUCENT **ORDER NUMBER: Z1.194**



DIMENSION (L x W):
330 x 250 mm

- | For X-ray of the head or the whole spinal column in inverse position
- | Mounting at foot end
- | Detachable padding
- | Angle adjustment by cogging

Suitable for:
» VARIUS

Max. head weight: 20 kg
Max. patient weight for inverse set-up 200 kg

MICROSURGERY RING, WITH INTEGRATED MOUNTING BRACKET

ORDER NUMBER: Z1.240




DIMENSIONS:
INTERIOR Ø: **325 mm**
EXTERIOR Ø: **440 mm**
PADDING WIDTH: **40 mm**

- | For precise operating and operating without tiring
- | Ideal for microscope surgery
- | Ball-and-socket hinge, height adjustable
- | Soft padding

Suitable for:
» FOCUS » VARIUS
» PRIMUS » GENIUS

MICROSURGERY WRIST REST (GENIUS EYE)

ORDER NUMBER: Z1.257



DIMENSIONS:

INTERIOR Ø: **292 mm**

EXTERIOR Ø: **372 mm**


- ▮ For precise operating and operating without tiring
- ▮ Ideal for microscope surgery
- ▮ Ball-and-socket hinge, height adjustable

Suitable for:

- » GENIUS EYE
- » GENIUS EYE PRO

LOWER ARM SUPPORT FOR OPHTHALMOLOGY

ORDER NUMBER: X.000013



- ▮ Ergonomic lower armrest for the surgeon during eye surgery in 12 o'clock position
- ▮ For operating without tiring during ophthalmological procedures
- ▮ Soft padding

Suitable for:

- » GENIUS EYE
- » GENIUS EYE PRO

TONGUE HOLDER, EXTENDED, WITH FLAT PLATE RAIL

ORDER NUMBER: Z1.061



DIMENSIONS (W x H):
460 x 390 – 670 mm

OPTIONAL: **TONGUE HOLDER, SHORT, WITH RAISED PLATE RIM (AS SHOWN)**

ORDER NUMBER: Z1.061

- ▮ For ENT
- ▮ Large individual range by way of two ball-and-socket hinges
- ▮ Chain for tongue and retractor holder
- ▮ instrument tray
- ▮ Height-adjustable

Suitable for:

- » FOCUS » VARIUS
- » PRIMUS » GENIUS

ARMREST OP WITH FASTENING BELT



ORDER NUMBER: Z1.001

PADDING DIMENSIONS (L x W x H):
490 x 140 x 40 mm

RESTRAINING BELT:
650 x 50 mm

- ▮ For ergonomic arm positioning
- ▮ Eccentric lever operates two ball-and-socket hinges, for fast, seamless adjustment
- ▮ Also for prone position
- ▮ Integrated impact protection
- ▮ Detachable padding, Silicone restraining belt

Suitable for:
» all BRUMABA Operating Tables

ARM SUPPORT LATERAL WITH FASTENING BELT



ORDER NUMBER: Z1.002

PADDING DIMENSIONS (L x W x H):
240 x 140 x 40 mm

- ▮ For slender arm apposition
- ▮ Folds down by way of swivel joint
- ▮ Detachable padding, Silicone restraining belt

Suitable for:
» all BRUMABA Operating Tables

PATIENT HANDLE, HINGED



ORDER NUMBER: Z1.033

PATIENT HANDLE, HINGED - GENIUS
ORDER NUMBER: Z1.281

PATIENT HANDLE, HINGED - VENUS
ORDER NUMBER: Z1.297

- ▮ For keeping the patient in the prone position

Suitable for:
» FOCUS » VARIUS
» PRIMUS

HAND SURGERY TABLE WITH MOUNTING BRACKET SPECIAL "XL"

ORDER NUMBER: Z1.037



PADDING DIMENSION (L x W):
660 x 400 mm

- ! For ergonomic positioning of the patient's arm
- ! Ball-and-socket hinge
- ! Radiolucent
- ! Slim design for maximum legroom

Suitable for:

- » FOCUS » VARIUS » VENUS
- » PRIMUS » GENIUS

HAND TABLE - TELESCOPIC SUPPORT, SINGLE

ORDER NUMBER: Z1.089



HEIGHT ADJUSTABLE FROM:
570 - 1030 mm

- ! For greater stability requirements at the hand table **Order number: Z1.037**
- ! Height-adjustable

Suitable for:

- » FOCUS » VARIUS » VENUS
- » PRIMUS » GENIUS

HAND SURGERY TABLE WITH DOUBLE TELESCOPIC SUPPORT, WITH MOUNTING BRACKET SPECIAL "XL"

ORDER NUMBER: Z1.036



PADDING DIMENSION (L x W):
800 x 500 mm

HEIGHT ADJUSTABLE FROM:
690 - 1110 mm

- ! For hand surgery with greater space requirements
- ! Radiolucent
- ! Ball-and-socket hinge
- ! Two height-adjustable telescopic supports for extra steadiness

Suitable for:

- » FOCUS » VARIUS » VENUS
- » PRIMUS » GENIUS

BACK SECTION CUSHION WITH ARM SUPPORT, VISCO

ORDER NUMBER: Z1.259



PADDING WIDTH:
720 mm


- For simple arm positioning

Suitable for:

- » FOCUS » VARIUS
- » PRIMUS » GENIUS

HORSESHOE-SHAPED HEADREST FOR SHOULDER ARTHROSCOPY WITH HEAD RESTRAINING BELT

ORDER NUMBER: Z1.154



PADDING DIMENSIONS (L x W x H):
200 x 210 x 80 mm

CLAMPING PIECE HEIGHT: **110 mm**

- For laterally offset positioning of the head
- For free access to the left and right shoulder area
- In combination with the head cradle, the horsehoe-shaped headrest adapts automatically during adjustment
- Ball-and-socket hinge, height adjustable, silicone restraining belt

Suitable for:

- » FOCUS » VARIUS
- » PRIMUS

BACK SECTION EXTENSION FOR SHOULDER ARTHROSCOPY

ORDER NUMBER: Z1.021



PADDING DIMENSIONS (L x W x H):
270 x 150 x 60 mm

- For shoulder operations left and right
- Easy assembly

Suitable for:

- » FOCUS » VARIUS
- » PRIMUS

SIDE GUARD, HINGED

ORDER NUMBER: Z1.187



DIMENSIONS (W x H):
350 x 310 mm

- ▮ For additional safety during transport of the patient
- ▮ Folds down by way of swivel joint

Suitable for:

» all BRUMABA Operating Tables

LATERAL POSITIONING SUPPORT, DEEP

ORDER NUMBER: Z1.110



PADDING DIMENSIONS (L x W x H):
180 x 300 x 50 mm

CLAMPING PIECE HEIGHT: **120 mm**

OPTIONAL:

LATERAL POSITIONING SUPPORT
ORDER NUMBER: Z1.054

CLAMPING PIECE HEIGHT: **65 mm**

- ▮ For stable lateral positioning of the patient
- ▮ Ball-and-socket hinge for seamless adjustment

Suitable for:

» FOCUS » VARIUS
» PRIMUS » GENIUS

LATERAL POSITIONING SUPPORT, DORSAL (BACK)

ORDER NUMBER: Z1.152



PADDING DIMENSIONS (L x W x H):
240 x 140 x 40 mm

HEIGHT-ADJUSTABILITY **120 mm**

- ▮ For slender support of the patient in lateral position
- ▮ Two joints for individual adaptation
- ▮ Height-adjustable

Suitable for:

» FOCUS » VARIUS
» PRIMUS » GENIUS



BODY BELT SPLITTED, SILICONE WITH 2 MOUNTING GLIDERS

ORDER NUMBER: Z1.241



DIMENSION (L x W):
1475 x 80 mm

- For securing the patient
- Hygienic silicon material
- Reduced pressure load

Suitable for:
» all BRUMABA Operating Tables

BODY BELT, SILICONE WITH 2 MOUNTING GLIDERS

**ORDER NUMBER: Z1.112
E.001061 | GENIUS EYE**



DIMENSION (L x W):
1440 x 80 mm

- For securing the patient
- Hygienic silicon material
- Reduced pressure load

Suitable for:
» all BRUMABA Operating Tables

LEG HOLDER OP

ORDER NUMBER: Z1.087



PADDING DIMENSIONS (L x W x H):

220 x 330 X 90 mm

RESTRAINING BELT DIMENSIONS (L x B):

650 x 50 mm

HEIGHT-ADJUSTABILITY: **130 mm**

- ▮ For ergonomic leg positioning
- ▮ Three-dimensional mobility
- ▮ Easy to handle by just one person
- ▮ Rocker function for automatic adaptation of the leg band against positional damage
- ▮ Ball-and-socket joints, height adjustable, silicon restraining belt

Suitable for:

- » FOCUS » VARIUS » VENUS
- » PRIMUS » GENIUS

MOUNTING ADAPTER WITH 90° DIRECTION CHANGE

ORDER NUMBER: X.00089



- ▮ For deflection of the surgical leg holder
- Order number: Z1.087**
for side positioning

Suitable for:

- » FOCUS » GENIUS
- » PRIMUS » VENUS
- » VARIUS

LEG HOLDER WITH 2 BALL-AND-SOCKET JOINTS

ORDER NUMBER: X.000126



PADDING DIMENSIONS (L x W x H):


220 x 330 X 90 mm

HEIGHT-ADJUSTABILITY **170 mm**

- ▮ For ergonomic leg positioning
- ▮ Three-dimensional mobility
- ▮ Easy to handle by just one person
- ▮ Suitable for side positioning
- ▮ Two ball-and-socket joints
- ▮ Height-adjustable
- ▮ Silicone restraining belt

Suitable for:

- » FOCUS » VARIUS » VENUS
- » PRIMUS » GENIUS



JOINT SUPPORT

ORDER NUMBER: Z1.028



PADDING DIMENSIONS (Ø x L):
90 x 230 mm

HEIGHT ADJUSTABLE FROM:
380 - 480 mm

- ! For surgery on the foot or shin
- ! Ball-and-socket joint
- ! Height-adjustable

Suitable for:

- » FOCUS » VARIUS » VENUS
- » PRIMUS » GENIUS

ARTHROSCOPIC LEG SUPPORT

ORDER NUMBER: Z1.016



MAX. Ø THIGH: **270 mm**


INCL. MOUNTING BRACKET XL
ORDER NUMBER: Z1.088

INCL. MOUNTING BRACKET
ORDER NUMBER: Z1.047

- ! for knee arthroscopic surgery
- ! Suitable for left and right knee
- ! Simple thigh adaptation
- ! Also for cruciate ligament surgery
- ! Can be rotated 360°
- ! Height-adjustable

Suitable for:

- » VARIUS



**PADDED INSERT FOR ARTHROSCOPIC SURGERY
LEG SUPPORT**

ORDER NUMBER: Z1.155




PADDING THICKNESS: 25 mm

■ For knee arthroscopic surgery
without tourniquet

Suitable for:
» VARIUS

STANDING PLATE (1 PAIR)

ORDER NUMBER: Z1.057



PADDING DIMENSIONS (B x L):
180 x 300 mm

■ For procedures in the reverse
Trendelenburg position, e.g. phlebology

Suitable for:
» BRUMABA Operating Tables "VARIUS"
and "VENUS"

FOR "PRIMUS"
ORDER NUMBER: x.000359

POSITIONING PADDINGS

LORDOSIS SUPPORT PADDING

ORDER NUMBER: Z1.140




PADDING DIMENSIONS (L x W x H):
400 x 210 x 60 mm

- Soft material and hygienic, seamless cover
- For the support of the lordosis

Suitable for:
» all BRUMABA Operating Tables

HALF ROLL PADDING, SMALL

ORDER NUMBER: Z1.100



PADDING DIMENSIONS (L x W x H):
400 x 140 x 70 mm

- Soft material and hygienic, seamless cover
- Comfortable shape for various purposes

Suitable for:
» all BRUMABA Operating Tables

HALF ROLL PADDING, LARGE

ORDER NUMBER: Z1.101



PADDING DIMENSIONS (L x W x H):
555 x 210 x 120 mm

- Soft material and hygienic, seamless cover
- Various purposes

Suitable for:
» all BRUMABA Operating Tables

POSITIONING PADDINGS

STOMACH POSITIONING PADDING, TRAPEZOIDAL (1 PIECE)

ORDER NUMBER: Z1.098



PADDING DIMENSION (L x W):
220 x 300 mm

MAX. DEPTH: **110 mm**

- Soft material and hygienic, seamless cover
- For ergonomic stomach positioning to relieve the abdominal area

Suitable for:

» all BRUMABA Operating Tables

MATTRESS, VISCOELASTIC

ORDER NUMBER: Z1.171



PADDING DIMENSIONS (L x W x H):
1660 x 570 x 40 mm

- For long operations to prevent bedsores and for extra comfort
- Anti-static

Suitable for:

» all BRUMABA Operating Tables

MATTRESS, VISCOELASTIC, SHORT, WITH GYNECOLOGICAL CUTOUT

ORDER NUMBER: Z1.213



PADDING DIMENSIONS (L x W x H):
950 x 550 mm

- For long operations to prevent bedsores and for extra comfort
- Anti-static

Suitable for:

» all BRUMABA Operating Tables



Mail




Internet



Index



TOC



BACK SECTION EXTENSION, NARROW

**ORDER NUMBER: Z1.006 | OP
Z1.183 | GENIUS**



PADDING DIMENSION (H x W):
100 x 420 mm

- For correct positioning of the shoulders in large patients
- Included as standard with FOCUS, PRIMUS, VARIUS and GENIUS

Suitable for:

- » FOCUS » VARIUS
- » PRIMUS » GENIUS

BACK SECTION EXTENSION, WIDE

**ORDER NUMBER: Z1.005 | OP
Z1.184 | GENIUS**



PADDING DIMENSION (H x W):
150 x 420 mm

- For correct positioning of the shoulders in large patients

Suitable for:

- » FOCUS » VARIUS
- » PRIMUS » GENIUS

BACK SECTION EXTENSION, EXTRA WIDE

ORDER NUMBER: Z1.212



PADDING DIMENSION (H x W):
150 x 550 mm

- ▮ For correct positioning of the shoulders in large patients
- ▮ Extra wide

Suitable for:

- » FOCUS » VARIUS
- » PRIMUS » GENIUS

BACK SECTION EXTENSION (GENIUS EYE / GENIUS EYE PRO)

ORDER NUMBER: Z1.279



PADDING DIMENSION (H x W):
150 x 530 mm

- ▮ For correct positioning of the shoulders in large patients

Suitable for:

- » GENIUS EYE
- » GENIUS EYE PRO

FOOT SECTION EXTENSION, WIDE

**ORDER NUMBER: Z1.268 | OP
Z1.286 | GENIUS***



PADDING DIMENSION (H x W):
150 x 420 mm

- ▮ For particularly large patients
- ▮ With edge protection plate
- ▮ Additional mounting bracket required


Suitable for:

- » PRIMUS » GENIUS

*INCL. MOUNTING ADAPTER

WITH OUT. MOUNTING ADAPTER

ORDER NUMBER: Z1.264 | GENIUS



**FOOT SECTION EXTENSION / SIDE EXTENSION
NARROW OP**

**ORDER NUMBER: Z1.321 | OP*
Z1.291 | GENIUS***



DIMENSIONS (H x W):
100 x 420 mm

OPTIONAL:
WITH MOUNTING RAIL
ORDER NUMBER: M.000644

! For particularly large patients or
widening of the lying surface

Suitable for:
» all BRUMABA Operating Tables

*INCL. MOUNTING ADAPTER
WITH OUT. MOUNTING ADAPTER

**ORDER NUMBER Z1.267 | OP
Z1.263 | GENIUS**

SIDE EXTENSION, RADIOLUCENT

ORDER NUMBER: Z1.055



DIMENSION (L x W x H):
410 x 225 X 55 mm

! For evacuation of the extremities during
operations and X-rays

! Additional mounting bracket required

Suitable for:
» FOCUS » VARIUS » VENUS
» PRIMUS » GENIUS

LOWER LEG PLATE, COMBINED (VARIUS)

ORDER NUMBER: Z1.170



PADDING DIMENSIONS | VISKO
(L x W x H):
500 x 555 x 60 mm

- For use during inverse positioning e.g. with external mounting
Order number: Z1.194 for a homogeneous X-ray image
- Connecting plate for the two-part foot section

Suitable for:
» VARIUS

LOWER LEG PLATE, COMBINED, EXTENDED (VARIUS)

ORDER NUMBER: Z1.169



PADDING DIMENSIONS | VISKO
(L x W x H):
800 x 555 x 60 mm

- For use during inverse positioning
- Extension of the X-ray range
- For a homogeneous X-ray image

Suitable for:
» VARIUS

LOWER LEG PLATES, EXTENDED, 1 PAIR (VARIUS)

ORDER NUMBER: Z1.020




PADDING DIMENSIONS (L x W x H):
800 x 270 x 50 mm

- Extension of the leg plates e.g. after lithotomy position or knee surgery

Suitable for:
» VARIUS

MOTOR MOUNTING BRACKET / INSTRUMENT TABLES



TRAY FOR MOTOR, SWIFELING, ON COLUMN

ORDER NUMBER: X.000040



PLATE SIZE: 400 x 360 mm

Suitable for:
» all BRUMABA Operating Tables

INSTRUMENT TABLE, MOBILE

ORDER NUMBER: Z1.172




PLATE SIZE: 600 x 400 mm

OPTIONAL:
ORDER NUMBER: Z1.282
PLATE SIZE: 700 x 500 mm

HEIGHT ADJUSTABLE FROM:
860 - 1100 mm

- ! For the storage of instruments within easy reach
- ! Mobile with hygienic, smooth-running double castors
- ! One-handed height adjustment via gas spring
- ! Brake against rolling away

INSTRUMENT TABLES



INSTRUMENT TABLE ON COLUMN

ORDER NUMBER: Z1.275



PLATE SIZE: 600 x 400 mm

OPTIONAL:

ORDER NUMBER: Z1.278

PLATE SIZE: 700 x 500 mm

HEIGHT-ADJUSTABLE ABOVE
TABLE LEVEL: 300 – 540 mm

- ! For the storage of instruments within ergonomic reach above the patient
- ! One-handed adjustment via gas spring
- ! Plate 360° rotatable, arm 180° tiltable
- ! Left and right mounting possible

Suitable for:

» all BRUMABA Operating Tables

TRAY FOR INSTRUMENT TABLE, ALUMINIUM

ORDER NUMBER: Z1.245




PLATE SIZE: 565 x 365 mm

- ! For efficient surgery organization
- ! Removable plate for the preparation of your instruments
- ! For instrument table
Order number: Z1.172 / Z1.275

Suitable for:

- » FOCUS
- » PRIMUS
- » VARIUS
- » GENIUS
- » GENIUS EYE
- » GENIUS EYE PRO

ANESTHESIA / INFUSION STANDS



ANESTHETIC FRAME BENDABLE, FOLDABLE

ORDER NUMBER: Z1.273



LENGTH: 960 mm

- ▮ For the supply of oxygen to the patient's head
- ▮ External hygienic tube guide on clips for 6 mm and 18 mm tubes (recommendation)
- ▮ Spacer for surgical drape
- ▮ Flexible gooseneck for individual adjustment
- ▮ Folds down by way of swivel joint

Suitable for:

» all BRUMABA Operating Tables

ANESTHETIC FRAME, BENDABLE - WITH OXYGEN SUPPLY

ORDER NUMBER: Z1.248



LENGTH: 1010 mm

- ▮ For the supply of oxygen to the patient's head
- ▮ Oxygen tube connection at gooseneck with outlet at head
- ▮ Spacer for surgical drape
- ▮ Flexible gooseneck for individual adjustment

Suitable for:

» all BRUMABA Operating Tables

ANESTHESIA / INFUSION STANDS

ANESTHESIOLOGY TUBE HOLDER

ORDER NUMBER: Z1.049



LENGTH: 300 mm + CROWN

- For holding all standard tube diameters (similar to the Ulmer wheel)
- Flexible gooseneck for individual adjustment,

Suitable for:

» all BRUMABA Operating Tables

ANESTHESIOLOGY TUBE HOLDER WITH REMOTE CONTROL HOLDER

ORDER NUMBER: X.000322




- For anesthesia tube holder

Suitable for:

» all BRUMABA Operating Tables

ANESTHESIA / INFUSION STANDS



ANESTHETIC FRAME, UNIVERSAL

ORDER NUMBER: Z1.015



WIDTH: **680 - 1130 mm**
HEIGHT-ADJUSTABLE: **270 mm**

- For all standard covering systems in surgery
- Length and height adjustable
- Ball-and-socket hinge
- Hook for infusion bottle

Suitable for:
» all BRUMABA Operating Tables

INFUSION STAND WITH BALL-AND-SOCKET JOINT

ORDER NUMBER: Z1.038



HEIGHT: **1110 - 1680 mm**
WITH 4 HOOKS
MAX. PERMISSIBLE LOAD:
2 KG PER HOOK

- Mounting on rail system
- One-handed height adjustment
- Ball-and-socket hinge
- 4 bottle holders

Suitable for:
» all BRUMABA Operating Tables

INFUSION STAND ON THE COLUMN, HINGED

ORDER NUMBER: Z1.277



HEIGHT: 1110 - 1680 mm

WITH 4 HOOKS

MAX. PERMISSIBLE LOAD:
2 KG PER HOOK

- Mounting on the column
- Independent of the angle of inclination setting of the patient surface
- One-handed height adjustment
- Tilttable
- 4 bottle holders

Suitable for:

» all BRUMABA Operating Tables

INFUSION STAND - SWIVEL ARM, EXTENDED (360 mm)

ORDER NUMBER: X.000214



OPTIONAL TO **ORDER**
NUMBER: Z1.181

LENGTH: 360 mm

- For a larger tilting range of the infusion stand
ORDER NUMBER: Z1.181

Suitable for:

» all BRUMABA Operating Tables

INFUSION STAND MOUNTED ON INSTRUMENT TABLE

ORDER NUMBER: Z1.280




- Hinged on instrument table on column
ORDER NUMBER: Z1.275 / Z1.278

- One-handed height adjustment
- 2 bottle holders

Suitable for:

» all BRUMABA Operating Tables

CONTROLS



HAND REMOTE CONTROL TYPE IV

ORDER NUMBER: Z1.118



- Included as standard with delivery of the operating table
- Magnetic back
- Carry out all individual settings
- Retrieval of all 8 memory positions

Suitable for:

- » FOCUS
- » PRIMUS
- » VARIUS
- » GENIUS
- » GENIUS EYE
- » GENIUS EYE PRO

REMOTE CONTROL HOLDER ON INFUSION STAND

ORDER NUMBER: Z1.243



- For attaching the hand-held remote control to the infusion stand

Suitable for:

- » all BRUMABA Operating Tables

REMOTE CONTROL HOLDER, FLEXIBLE

ORDER NUMBER: Z1.149



GOOSENECK LENGTH : 300 mm


- ▮ Flexible gooseneck

Suitable for:

- » all BRUMABA Operating Tables

REMOTE CONTROL HOLDER, SHORT

ORDER NUMBER: Z1.026



- ▮ Mounting on rail system underneath patient surface

Suitable for:

- » all BRUMABA Operating Tables

MAGNETIC PLATE FOR REMOTE CONTROL

ORDER NUMBER: Z1.045




- ▮ Included as standard with delivery of every operating table

- ▮ Attachment by way of rear adhesive surface

Suitable for:

- » all BRUMABA Operating Tables

CONTROLS



FOOT CONTROL, TYPE IV

ORDER NUMBER: Z1.185.01
ORDER NUMBER: Z1.303 | VENUS



WEIGHT: 4.25 kg

- ▮ Carry out all individual settings
- ▮ Retrieval of all 8 memory positions
- ▮ Lifting bar for lifting or moving the foot control with the foot
- ▮ Can be operated in sterile condition

Suitable for:

- » FOCUS » GENIUS
- » PRIMUS » GENIUS EYE
- » VARIUS » GENIUS EYE PRO

FOOT CONTROL, TYPE IV WITHOUT MEMORY

ORDER NUMBER: Z1.271.01



WEIGHT: 4.25 kg

- ▮ Carry out all individual settings
- ▮ Lifting bar for lifting or moving the foot control with the foot
- ▮ No retrievable memory positions

Suitable for:

- » all BRUMABA Operating Tables

PLATFORM FOR FOOT SWITCH TYPE IV - 4-EASY DRIVE

ORDER NUMBER: Z1.274



- ! For positioning the foot control on the undercarriage 4-EASY-DRIVE, 5-WHEEL-DRIVE

Suitable for:

- » all BRUMABA Operating Tables

FOOT CONTROL, TYPE IV FIXED ON UNDERCARRIAGE

ORDER NUMBER: Z1.225.01



WEIGHT: 4.25 kg

- ! Carry out all individual settings
- ! Retrieval of all 8 memory positions
- ! Foot control fixed on undercarriage 4-EASY-DRIVE, 5-WHEEL-DRIVE

Suitable for:

- » FOCUS » GENIUS
- » PRIMUS » GENIUS EYE
- » VARIUS » GENIUS EYE PRO

FOOT CONTROL, TYPE IV EYE, WITHOUT MEMORY

ORDER NUMBER: Z1.226.01




WEIGHT: 4.25 kg

- ! Adjustment of head section, back section and height possible
- ! Foot control fixed on undercarriage
- ! No retrievable memory positions

Suitable for:

- » all BRUMABA Operating Tables



MOUNTING BRACKET

ORDER NUMBER: Z1.004



FLANGE DIMENSION (L x B):
30 x 10 mm

- ▮ For fastening of BRUMABA accessories
- ▮ Safe 4-point clamping

Suitable for:
» all BRUMABA Operating Tables

MOUNTING BRACKET SPECIAL "XL"

ORDER NUMBER: Z1.088



FLANGE DIMENSION (L x B):
30 x 15 mm

- ▮ For fastening of BRUMABA accessories
- ▮ Safe 4-point clamping

Suitable for:
» all BRUMABA Operating Tables

MOUNTING GLIDER WITH CLAMPING

ORDER NUMBER: Z1.093




FLANGE DIMENSION (L x B):
30 x 10 mm

! For the permanent fastening of accessories

Suitable for:
» all BRUMABA Operating Tables

MOUNTING ADAPTER, 25 x 10 mm

ORDER NUMBER: Z1.047



DIMENSIONS (L x H x B):
200 x 25 x 10 mm

OPTIONAL: **MOUNTING ADAPTER, 30 x 10 mm (SPECIAL)**

ORDER NUMBER: Z1.047A

DIMENSIONS (L x H x B):
200 x 30 x 10 mm

! For fastening accessories of other manufacturers (standard rails)

Suitable for:
» all BRUMABA Operating Tables

MOUNTING ADAPTER, STRENGTHENED, 290 mm

ORDER NUMBER: X.000244



DIMENSIONS (L x H x B):
290 x 25 x 10 mm

! For fastening accessories of other manufacturers (standard rails)

Suitable for:
» all BRUMABA Operating Tables

MOUNTING RAIL, 250 mm HEAD SECTION (GENIUS)

ORDER NUMBER: Z1.151



! For head section left

Suitable for:
» GENIUS

MOUNTING RAIL, 340 mm FOOT SECTION (GENIUS)

ORDER NUMBER: Z1.079



! For foot section left or right

Suitable for:
» GENIUS

MOUNTING RAIL, 250MM MASSIV SEAT SECTION (GENIUS)

ORDER NUMBER: X.000087



! For seat section, to attach leg holders
ORDER number: Z1.087

Suitable for:
» GENIUS

**MOUNTING RAIL, 340 MM,
SOLID BACK AND FOOT SECTION (GENIUS)**

ORDER NUMBER: X.000036




! For back and foot section

Suitable for:
» GENIUS

**MOUNTING RAIL, 125 mm
HEAD SECTION (GENIUS EYE / EYE PRO)**

ORDER NUMBER: Z1.262




! E.g. for mounting the
microsurgery ring GENIUS EYE
ORDER NUMBER: Z1.257

Suitable for:
» GENIUS EYE
» GENIUS EYE PRO

MOUNTING RAILS SET (VENUS)

ORDER NUMBER: Z1.304



! For back, seat and foot section

Suitable for:
» VENUS

TRANSPORT HANDLES

TRANSPORT HANDLES, FOLDAWAY, HEAD END, GENIUS (1 PAIR)



ORDER NUMBER: Z1.210

- For the back-friendly moving of the operating table in chair position
- Foldaway for legroom

Suitable for:

- » GENIUS
- » GENIUS EYE | included as standard
- » GENIUS EYE PRO | included as standard

TRANSPORT HANDLES, FOLDAWAY, HEAD END, FOCUS, PRIMUS, VARIUS (1 PAIR)



ORDER NUMBER: Z1.211

- For the back-friendly moving of the operating table in chair position

Suitable for:

- » FOCUS » VARIUS
- » PRIMUS

TRANSPORT HANDLES, HINGED ON FOOT SIDE (1 PAIR)

ORDER NUMBER: Z1.164



- For the back-friendly moving of the operating table from the foot end


Suitable for:

- » all BRUMABA Operating Tables

MOUNTING CLAMP XL REQUIRED

ORDER NUMBER: Z1.088

HYGIENE / PRESERVATION / CARE



EASYCOVER (150 PIECES - DISPOSABLE, UNSTERILE)

ORDER NUMBER: Z1.138



DIMENSION (L x W):
2000 x 720 mm

PACKAGING UNIT: **150 PIECES**
COLOR: **WHITE**

- ! Watertight disposable surgery bed drape for maximum hygiene
- ! Head pocket prevents cover from slipping
- ! Paper coating for comfort while lying

Suitable for:

» all BRUMABA Operating Tables

PAPER TOWEL HOLDER

ORDER NUMBER: Z1.051



DIMENSIONS (W):
560 mm

USE OF PAPER TOWELS UP TO
MAX. WIDTH OF 560 mm

- ! For attaching paper towels under the foot section

Suitable for:

- » FOCUS » GENIUS
- » PRIMUS » GENIUS EYE
- » VARIUS » GENIUS EYE PRO

VENUS ORDER NUMBER: Z1.294

DIMENSIONS (W): **550mm**

HYGIENE / PRESERVATION / CARE

PROTECTIVE COVER FOR FOOT SECTION, VISCO



ORDER NUMBER: Z1.223

- For the mechanical protection of the foot section padding e.g. outdoor footwear
- Easy disassembly for cleaning


Suitable for:

- » FOCUS
- » PRIMUS

Optional:

**ORDER NUMBER: Z1.224 | GENIUS
Z1.299 | VENUS**

PROTECTIVE COVER FOR FOOT SECTION, VISCO, (VARIUS) (1 PAIR)



ORDER NUMBER: Z1.148A

- For the mechanical protection of the foot section padding e.g. outdoor footwear
- Easy disassembly for cleaning

Suitable for:

- » VARIUS

SERVICE-KIT



ORDER NUMBER: Z1.305

- Everything to care for and clean your BRUMABA operating table
- Matting sponge, sponge, microfiber cloth, lubricant, cold welding paste, upholstery care, stainless steel care

Suitable for:

- » all BRUMABA products

LAMP WITH LOUPE ON COLUMN, HINGED

ORDER NUMBER: Z1.217



LENS DIMENSION Ø: **160 mm**

LENS DIOPTRER **3.5**

COLOR RENDITION: **6500 K**

LED ENERGY EFFICIENCY: **A+**

Suitable for:

- » FOCUS » GENIUS
- » PRIMUS » GENIUS EYE PRO
- » VARIUS

Optional:

ORDER NUMBER: Z1.300 | VENUS

ACCESSORIES ORDER

WALL MOUNTING RAIL



ORDER NUMBER: Z1.162

LENGTH: 100 mm

- For fastening accessories to the wall
- 5 gliders included

Suitable for:

» all BRUMABA Operating Table accessories

GLIDERS FOR WALL MOUNTING /
ACCESSORIES TROLLEY:

ORDER NUMBER: Z1.161

WALL MOUNT FOR HAND SURGERY TABLE



ORDER NUMBER: Z1.190

- For space-saving storage of the hand table

Suitable for:

Hand surgery table Z1.037

ACCESSORIES COLUMN



ORDER NUMBER: Z1.117

HEIGHT: 1080 mm

Ø BASE: 620 mm

- For storage of up to 8 accessories
- Mobile

Suitable for:

» all BRUMABA Operating Table accessories

Optional accessory:

STORAGE PLATE **ORDER NUMBER: Z1.153**

ACCESSORIES COLUMN - STORAGE PLATE



ORDER NUMBER: Z1.153

Ø STORAGE PLATE: 620 mm

■ Additional storage option

for the accessory column

ORDER NUMBER: Z1.117

Suitable for:

» all BRUMABA Operating Table accessories

ACCESSORIES TROLLEY



ORDER NUMBER: Z1.060

DIMENSIONS (L x W x H):

730 x 550 x 1000 mm

RAIL LENGTH: **725 mm**

■ For the storage of accessories

■ Mobile

■ 12 gliders included

■ Extension for hand table and arthroscopic surgery leg support

Suitable for:

» all BRUMABA Operating Table accessories

GLIDERS FOR WALL MOUNTING /
ACCESSORIES TROLLEY:

ORDER NUMBER: Z1.161

TRANSPORT TROLLEY FOR BATTERIES



ORDER NUMBER: Z1.175

DIMENSIONS (L x W x H):

300 x 450 x 800 mm

■ For easy transportation of two battery packs

■ Holder for the charger

ORDER NUMBER: Z1.134

Suitable for:

» FOCUS » VARIUS

» PRIMUS

INDEX

A

Accessories column	46
Accessories column - storage plate	47
Accessories trolley	47
Anesthetic frame bendable, foldable	30
Anesthetic frame, universal	32
Anesthetic frame, bendable - with oxygen supply	30
Anesthesiology tube holder	31
Anesthesiology tube holder with remote control holder	31
Armrest OP with fastening belt	14
Arm support lateral with fastening belt	14
Arthroscopic leg support	20

B

Back section cushion with arm support, visco	16
Back section extension (GENIUS EYE / GENIUS EYE PRO)	25
Back section extension, extra wide	25
Back section extension for shoulder arthroscopy	16
Back section extension, narrow (GENIUS)	24
Back section extension, narrow (OP)	24
Back section extension, wide (GENIUS)	24
Back section extension, wide (OP)	24
Body belt, silicone with 2 mounting gliders	18
Body belt splitted, silicone with 2 mounting gliders	18

E

EasyCover (150 pieces - disposable, unsterile)	43
External mounting with head plate, radiolucent	12

F

Foot control, type IV	36
Foot control, type IV fixed on undercarriage	37
Foot control, type IV EYE, without memory	37
Foot control, type IV without memory	36
Foot section extension, wide (GENIUS)	25
Foot section extension, wide (OP)	25

G

Gliders for wall mounting rail / accessories trolley	46
---	----

H

Hair transplant support	10
Hand remote control type IV	34
Half roll padding, small	22
Half roll padding, large	22
Hand surgery table with mounting bracket special "XL"	15
Hand surgery table with double telescopic support, with mounting bracket special "XL"	15
Hand table telescopic support, single	15
Head plate, large, 400 x 300 mm	8
Head plate, small, 320 x 240 mm	8
Head plate "rubins pillow" with ball-and-socket joint	10
Headrest, detachable	11
Headrest, detachable, GENIUS (surcharge)	11
Horseshoe-shaped headrest, (included as standard)	8
Horseshoe-shaped headrest for shoulder arthroscopy with head restraining belt	16
Horseshoe-shaped headrest with long clamping piece, 110 mm	9
Horseshoe-shaped headrest, small	9

I

Infusion stand mounted on instrument table	31
Infusion stand on the column, hinged	33
Infusion stand - swivel arm, extended (360 mm)	33
Instrument table, mobile, 600 x 400 mm	28
Instrument table, mobile, 700 x 500 mm	28
Instrument table on the column, 600 x 400 mm	29
Instrument table on the column, 700 x 500 mm	29

J

Joint support	20
---------------	----

L

Lateral positioning support, deep	17
Lateral positioning support, dorsal (back)	17
Leg holder with 2 ball-and-socket joints	19
Leg holder, OP	19
Lordosis support padding	22
Lower arm support for ophthalmology	13
Lower leg plate, combined, extended (VARIUS)	27



Lower leg plate, combined (VARIUS) _____ 27
 Lower leg plates, extended, 1 pair (VARIUS) _____ 27

M

Magnetic plate for remote control _____ 35
 Magnifying lamp with swivel arm _____ 45
 Magnifying lamp with swivel arm (VENUS) _____ 45
 Mattress, viscoelastic _____ 23
 Mattress, viscoelastic, short, with gynecological cutout _____ 23
 Microsurgery ring, with integrated mounting bracket _____ 12
 Microsurgery wrist rest (GENIUS EYE) _____ 13
 Mounting adapter with 90° direction change _____ 19
 Mounting adapter, 25 x 10 mm _____ 39
 Mounting adapter, strengthened, 290 mm _____ 39
 Mounting bracket _____ 38
 Mounting bracket special "XL" _____ 38
 Mounting glider with clamping _____ 39
 Mounting rail, 250 mm head section (GENIUS) _____ 40
 Mounting rail, 340 mm foot section (GENIUS) _____ 40
 Mounting rail, 250 mm, massiv seat section (GENIUS) _____ 40
 Mounting rail, 340 mm, solid back and foot section (GENIUS) _____ 41
 Mounting rail, 125 mm head section (GENIUS EYE, -EYE PRO) _____ 41
 Mounting rail set (VENUS) _____ 41

N

Neck cushion, viscoelastic, black _____ 10

P

Padded insert for arthroscopic surgery leg support _____ 20
 Patient handle, hinged _____ 14
 Paper towel holder _____ 43
 Paper towel holder, 550 mm (VENUS) _____ 43
 Platform for foot switch type IV - 4-EASY DRIVE _____ 37
 Protective cover for foot section, visco _____ 44
 Protective cover for foot section, visco, (VARIUS) (1 pair) _____ 44
 Side extension, radiolucent _____ 26
 Side guard, hinged _____ 17
 Silicone head restraining belt for fastening on the headrest _____ 9
 Standing plate (1 pair) _____ 20
 Stomach positioning padding, trapezoid _____ 23

R

Radiolucent headrest, detachable (surcharge) _____ 11
 Remote control holder, flexible _____ 35
 Remote control holder on infusion stand _____ 34
 Remote control holder, short _____ 35

S

Service kit _____ 44
 Side widener / foot section extension, narrow (GENIUS) _____ 26
 Side extension / foot section extension, narrow (OP) _____ 26
 Side extension, narrow _____ 26
 Side extension, radiolucent _____ 26
 Side guard, hinged _____ 17
 Silicone head restraining belt for fastening on the headrest _____ 9
 Standing plate (1 pair) _____ 20
 Stomach positioning padding, trapezoid _____ 23

T

Tongue holder, extended, with flat plate rail _____ 13
 Transport handles, foldaway, head end (1 pair) (GENIUS) _____ 42
 Transport handles foldaway, head end, (1 pair) (FOCUS, PRIMUS, VARIUS) _____ 42
 Transport handles, hinged on foot side (1 pair) _____ 42
 Transport trolley for batteries _____ 47
 Tray for motor, swifeling, on column _____ 28
 Tray for instrument table, aluminium _____ 29

W

Wall mount for hand surgery table _____ 46
 Wall mounting rail _____ 46





Mail



Internet



Index



TOC



ADRESSE

BRUMABA GmbH & Co. KG
Bürgermeister-Graf-Ring 17
82538 Geretsried | Germany



TELEFON & FAX

Tel: +49 (0) 8171 / 2672 - 0
Fax: +49 (0) 8171 / 2672 - 10

ONLINE

Email: info@brumaba.de
Web: www.brumaba.de

© BRUMABA GmbH & Co. KG | 2020 | 0120-R-250.2