



# Diathermy

Endodiathermy ..... 36  
Exodiathermy ..... 37



## Endodiathermy

D8203

### Intraocular Diathermy Probe



Designed for retinal endodiathermy. 31 mm shaft tapered from 21GA to 23GA at very delicate tip. Plastic covered handle. Autoclavable. Lemosa type adapter. Overall length 132 mm.



E7920

### Bipolar Eraser, 25GA, tapered



Designed for retinal endodiathermy. 25GA shaft. Conical tip. Plastic covered handle. Autoclavable. Lemosa type connector.



E7924

### Bipolar Eraser, 25GA, 2pin tapered



Designed for retinal endodiathermy. 25GA shaft. Conical tip. 2-pin connector. Sterile single-use. 12/box.



E7925

### Bipolar Eraser, 25GA, 2pin blunt



Designed for retinal endodiathermy. 25GA shaft. 45° beveled tip. 2-pin connector. Sterile single-use. 12/box.



E7926

### Bipolar Eraser, 25GA, lemo tapered



Designed for retinal endodiathermy. 25GA shaft. Conical tip. Lemosa type connector. Sterile single-use. 12/box.



E7927

### Bipolar Eraser, 25GA, lemo blunt



Designed for retinal endodiathermy. 25GA shaft. 45° beveled tip. Lemosa type connector. Sterile single-use. 12/box.



E7929

### Bipolar Eraser, 20GA, lemo tapered



Designed for retinal endodiathermy. 20GA shaft. Conical tip. Lemosa type connector. Sterile single-use. 12/box.



## Endodiathermy

E7930

Bipolar Eraser, 23GA, lemo tapered



Designed for retinal endodiathermy. 23GA shaft. Conical tip. Lemos type connector. Sterile single-use. 12/box.



E7931

Bipolar Eraser, 23GA, 2pin tapered



Designed for retinal endodiathermy. 23GA shaft. Conical tip. 2-pin connector. Sterile single-use. 12/box.



E7932

Bipolar Eraser, 27GA, lemo tapered



Designed for retinal endodiathermy. 27GA shaft. Conical tip. Lemos type connector. Sterile single-use. 12/box.



E7933

Bipolar Eraser, 27GA, 2pin tapered



Designed for retinal endodiathermy. 27GA shaft. Conical tip. 2-pin connector. Sterile single-use. 12/box.



## Exodiathermy

D8200

Straight Bipolar Forceps



0.5 mm tip. Autoclavable. Lemos type adapter. Overall length 90 mm.



D8201

Bipolar Eraser



For scleral haemostasis. 28.5 mm shaft with 45° beveled tip. Plastic covered handle. Autoclavable. Lemos type adapter. Overall length 131 mm.



## Exodiathermy

E7928

### Bipolar Scleral Eraser




Designed for scleral haemostasis. 45° beveled tip. Lemossa type connector. Sterile single-use. 12/box.



S20501

### Jeweler Type Bipolar Forceps



Straight shafts with micro tips and dull finish. Smooth handle with dull finish. Overall length 119 mm.



S205010A

### McPherson Angled Bipolar Forceps




Straight shafts with 4 mm jaws. Serrated handle with dull finish. Jaws angled 45°. Overall length 91 mm.



S205010S

### McPherson Straight Bipolar Forceps



Straight shafts with 4 mm jaws. Serrated handle with dull finish. Jaws straight. Overall length 94 mm.



S205014A

### Shepard Bipolar Forceps Angled



Angled shafts with 4 mm platform. Wide serrated handle with dull finish. Platform angled 45°. Overall length 115 mm.

