



# FLUO 900

Safe, simple and efficient

100 STERILE PAPER STRIPS

FLUO 900

## Safe, simple and efficient

### Safe

The individually wrapped sterile paper strip ensures safe application of the fluorescein.

### Simple

Just moisten the strip with sterile ophthalmic solution and touch conjunctiva or fornix with moistened tip to dispense the dye.

### Efficient

With 1 mg Fluorescein Sodium I.P. per strip, one strip contains sufficient dye to obtain optimal staining results.



Distributed by:  
HAAG-STREIT AG  
Gartenstadtstrasse 10  
3098 Koeniz  
Switzerland  
Phone +41 31 978 01 11  
Fax +41 31 978 02 82  
info@haag-streit.com  
www.haag-streit.com

Manufactured by:  
Entod Research Cell UK Ltd.  
North London Business Park  
Oakleigh Road South, New  
Southgate  
London N11 1GN, UK