

Microsurgical Accesories

A variety of products to complement ophthalmic microsurgery and minor office procedures.


Fluid Absorption and Administration




- 9001 Eye Spear**
Fast wicking action for quick absorption



- 9003 Absorbent Stick - Sterile**
2 ³/₈ in (6.6cm) overall length ³/₁₆ in (5mm) diameter, 45° angled ends. Soft, non-woven material is rigid enough to withstand gentle pressure. Both ends may be used for absorption



- 9004 Absorbent Stick - Non-Sterile**
2 ³/₈ in (6.6cm) overall length ³/₁₆ in (5mm) diameter, 45° angled ends. Soft, non-woven material is rigid enough to withstand gentle pressure. Both ends may be used for absorption



- 9006 Surgical Sponge - Gills**
2 x 8mm - Lint free PVA, useful for applying or removing fluid in a precise location




- 9101 Instrument Wipe**
3 x 3 in (7.5 x 7.5cm) Lint free PVA for cleaning delicate instruments. Soft when moistened and will not drag on blades or tips



- 9301 .20 Micron Filter**
Filters out particles and microbiological contamination .22 micron and larger. Male/female luer lock fittings




- 9303 Male to Male Connector**
Overall length 29mm
Luer slip fittings connect instrument to tubing
Especially useful for bimanual I and A to connect cannula to I/A tubing




- 9701 Incise Film**
Overall size: 5 x 9 in (12.7 x 22.8cm)
Adhesive area: 5 x 7 ¹/₄ in (12.7 x 18.4cm)
Small incise covers lashes, lids and brow



- 9802 Eye Fluid Wick**
Highly absorbent lint free PVA wick
Removes fluid and directs it away from surgical area



- 9803 Surgical Absorbent Kit**
1 - 9101 Instrument Wipe
1 - 9202 Corneal Light Shield
1 - 9802 Eye Fluid Wick
Popular microsurgical accessories combined into one package for convenience and cost savings




- 7002 Silicone Tubing**
Length: 14 in (35cm)
ID: .012 in (.30mm)
OD: .025 in (.60mm)

Microsurgical Accessories

A variety of products to complement ophthalmic microsurgery and minor office procedures.

Eye Protection Protects the eye during or after surgery or trauma.



9401 Clear Eye Shield - Non-sterile
Smooth edges need no protective cover. Universal for left or right eye




9402 Clear Eye Shield - Sterile
Smooth edges need no protective cover. Universal for left or right eye



9202 Corneal Light Shield
8mm diameter - Lint free PVA
Keeps cornea moist. Protects retina from operating microscope light rays




9203 Corneal Light Shield
6mm diameter - Lint free PVA
Keeps cornea moist. Protects retina from operating microscope light rays



9205 Corneal Light Shield
10mm diameter - Lint free PVA
Keeps cornea moist. Protects retina from operating microscope light rays

Surgical Glides Useful for placement of intraocular lenses into the posterior chamber. Also protects iris from prolapse or trauma during microsurgical procedures.



9201 Surgical Glide
25mm long x 3mm wide
1.5mm radius; 6 mil thick



9204 IOL Glide - Sheets
35mm long x 5mm wide
2.5mm radius; 6 mil thick

