



Probes

Lacrimal, DCR 152



Lacrimal, DCR

51011

Wilder Lacrimal Dilator

BAUSCH + LOMB
Instruments



Fine - width at tip 0.2 mm. Tapered shaft to 3.5 mm.
Shaft length 25 mm. Overall length 101 mm.



51012

Wilder Lacrimal Dilator

BAUSCH + LOMB
Instruments



Medium - width at tip 0.4 mm. Tapered shaft to 3.5 mm.
Shaft length 25 mm. Overall length 101 mm.



51013

Wilder Lacrimal Dilator

BAUSCH + LOMB
Instruments



Strong - width at tip 0.7 mm. Tapered shaft to 3.5 mm.
Shaft length 25 mm. Overall length 101 mm.



51100

Lacrimal Dilator, Model Giessen

BAUSCH + LOMB
Instruments



Overall length 85 mm.



51130

Lacrimal Dilator

BAUSCH + LOMB
Instruments



Fine. Overall length 75 mm.



51171

Pigtail Probe











BAUSCH + LOMB
Instruments



With 0.5 mm suture eye. Overall length 145 mm.



Lacrimal, DCR


51171Y1	Pigtail Probe	BAUSCH + LOMB Instruments	
<p>0.5 mm suture eye. Child size. Overall length 140 mm.</p>			
E0421	Harms Trabeculotome Set	BAUSCH + LOMB Instruments	
<p>9 mm long pointed tips with 3 mm spread. Pair consists of left and right. Round knurled handle. Polished finish. Overall length 46 mm.</p>			
E4198	Bowman Stainless Probe - 0000 and 000	BAUSCH + LOMB Instruments	
<p>Sizes 0000 and 000. Malleable stainless steel. Polished finish. Overall length 141 mm.</p>			
E4251	Worst Pigtail Probe	BAUSCH + LOMB Instruments	
<p>9 mm diameter pigtail curved blunt probes. Includes pre-cut polyethylene tubes (E4251 B). Overall length 149 mm.</p>			
E4251C	Pigtail Probe	BAUSCH + LOMB Instruments	
<p>8 mm diameter pigtail curved probes with suture eyes just back of blunt ends. Overall length 143 mm.</p>			

Lacrimal, DCR

E4340

Castroviejo Lacrimal Dilator

BAUSCH + LOMB
Instruments



One end with needle point and other end with medium taper. Overall length 140 mm.



E4366

Ruedemann Lacrimal Dilator

BAUSCH + LOMB
Instruments



12 mm short taper. Overall length 71 mm.

