



Scleral Depressors



Scleral Depressors

56040A

Meyer-Schwickerath Localizer

Overall length 123 mm.



56201A

Laqua Scleral Depressor


Double-ended. 4 mm teardrop-shaped end with 34 mm shaft and 2 mm trephine-like blunt end with 33 mm shaft. Overall length 136 mm.



56210A

Kreissig Localizer

Double-ended. 2.4 mm diameter conical tip. Slightly curved 32 mm shaft. 6.5 mm wide bar tip. Strongly curved 35 mm shaft. Overall length 135 mm.



E5100 10

Schepens Scleral Depressor Medium

Polished depressor with 8 mm wide tip. Dull finished thimble. Size 10, medium thimble. Overall length 43 mm.



E5107

Flynn Scleral Depressor


Lens loop style tip ideal for depressing pediatric sclera. Flat serrated handle. Polished finish. Overall length 80 mm.



E5108

Schocket Double Ended Scleral Depressor

Double-end depressor and marker. Rounded end designed for depressing sclera and curved marking end for reaching behind the globe. 6 mm wide bar and 2.5 mm teardrop shape tips. Overall length 134 mm.




E5109

O'Connor Scleral Depressor-Marker




3 mm teardrop shaped end and 1.5 mm Gass marker end. Flat, serrated handle. Polished finish. Overall length 133 mm.



E5112

Spaide Pediatric Scleral Depressor



Tip curved to fit the globe of an infant. Serrated surface grasps the conjunctiva to rotate the eye in any direction. The tip of the foot is used to depress the sclera. Flat serrated handle. Dull finish. Overall length 125 mm.



