

St
Stellaris

Ak
Akreos

In
Instruments

Vi
Viscoelastic

Su
Support

The MICS™ Platform

The Exclusive 1.8 mm Solution From Bausch & Lomb

Stellaris® MICS

The Stellaris Vision Enhancement System delivers the critical capabilities required for 1.8 mm MICS: superior fluidics control, efficient cutting dynamics and the ergonomic insights that make the system easy to use.

Akreos® MICS

The Akreos MICS Lens combines the vision benefits of an aspheric acrylic lens with the material and design characteristics that allow it to be inserted into a true 1.8 mm incision.

MICS Instruments

The MICS Platform includes precision performance instruments from STORZ® Ophthalmics for both C-MICS and B-MICS procedures.

MICS Viscoelastic

Amvisc® PLUS is the versatile viscoelastic that is ideal for every step of your 1.8 mm MICS procedure.

MICS Support

When you move to the MICS Platform you benefit from a comprehensive support program designed to assist cataract surgeons in developing and perfecting their MICS techniques.

St
Stellaris

Innovative Technology
For 1.8 mm Micro Incision Cataract Surgery



Stellaris® MICS™ Vision Enhancement System

Safe, Efficient, Easy to Use

Choose

StableChamber™ Fluidics allows customization based on your fluidics preference, choosing either flow or vacuum control.

BL11110
Vacuum Fluidics Module (VFM)



BL11120
Advanced Fluidics Module (AFM)



Foot Control

Stellaris® wireless Dual-Linear foot control




Options Available

BL2111
Stellaris® remote control



Note: Button positions correspond to icon positions on screen

BL2300
Stellaris® multi media center (MMC)



Video inlay





Video overlay

BL6320
DigiFlow™

The DigiFlow software upgrade enables control of automatic pressurization of the BSS bottle, minimizing IOP fluctuations during cataract surgery.
Please ask your Bausch & Lomb representative for details

DigiFlow™ requires use of:
D4600A
DigiFlow™ Tubing



BL3170
Stellaris® Attune™ ultrasonic handpiece



The six-crystal handpiece delivers optimized 28.5 kHz frequency for consistent power and increased stroke length.

BL3160
Stellaris® I/A handpiece



*Used with the BL31xx series of I/A tips.
Proprietary locking irrigating connector for use with the Stellaris™ System tubing kits.*

1.8 mm MICS Solution

"The conversion from standard phaco to Stellaris 1.8 mm MICS was a breeze. The safety and control with 1.8 MICS was there from the first case on... Because of the advancements with Stellaris 1.8 MICS, I exclusively use this technology today." - John D. Hunkeler, MD

Single Use Knives

6 per box **STERILE**

E7596

1.5-2.0 mm
trapezoidal



E7600

1.6-1.8 mm
trapezoidal



E7602

1.8-2.2-2.8 mm
uni-trap



E7515

15° sideport



Capsulorhexis Forceps

ET1088TS

Capsulorhexis
forceps tip



ET8190 H

Capsulorhexis forceps
handle titanium (non-rotating)



ET1202

Utrata Micro Capsulorhexis
Forceps



ET1203T

MICS Capsulorhexis Forceps
Tip with Seibel Ruler



ET1203H

MICS Capsulorhexis Forceps
Handle for E2013T




Phaco Packs Contents

6 per box **STERILE**

BL5113A

StableChamber™ Phaco Pack with MICS™ needle



1 VFM cassette assembly
and I/A set
1 Auxiliary drape
1 Screen drape
1 Remote control drape
2 Irrigation sleeves
1 Plug
1 Test chamber
1 Needle & wrench

BL5123A

Micro Incision AFS Phaco Pack with MICS™ needle



1 AFM cassette assembly
and I/A set
1 Auxiliary drape
1 Screen drape
1 Remote control drape
2 Irrigation sleeves
1 Test chamber
1 Needle & wrench

Phaco Packs	Stable Chamber	Phaco Needle C-MICS 1.8mm	
		RIS* 1.8 mm OD* 0.74-0.95 mm ID* 0.50-0.79 mm	
		Straight	Angled 30°
VFM	BL5113	•	•
	BL5113A	•	•
	BL5114	•	•
	BL5114A	•	•
AFM	BL5123	•	•
	BL5123A	•	•

* RIS - Recommended Incision Size - * OD - Outer Diameter - * ID - Inner Diameter

C-MICS Coaxial I/A Single Use

12 per box **STERILE**

85794ST

I/A handpiece curved
shaft tip, with sleeve
for 1.8 mm incision



85795ST

I/A handpiece
45° angled tip, with
sleeve for 1.8 mm incision



2.2 mm Small Incision Solution

Single Use Knives

E7601

2.0-2.2 mm trapezoidal



E7602

1.8-2.2-2.8 mm uni-trap



E7515

15° sideport



6 per box **STERILE**

Capsulorhexis Forceps

ET1202

Utrata Micro Capsulorhexis Forceps



ET2012

Devgan Micro Capsulorhexis Forceps




Phaco Packs Contents

6 per box **STERILE**

BL5112

Optimized Stability Vacuum Pack



- 1 VFM cassette assembly and I/A set
- 1 Auxiliary drape
- 1 Screen drape
- 1 Remote control drape
- 2 Irrigation sleeves
- 1 Plug
- 1 Test chamber
- 1 Needle wrench

BL5120

Premium AFS Phaco Pack



- 1 AFM cassette assembly and I/A set
- 1 Auxiliary drape
- 1 Screen drape
- 1 Remote control drape
- 2 Irrigation sleeves
- 1 Test chamber
- 1 Needle wrench

Phaco Needles

6 per box **STERILE**



Single use	MicroFlow® 2.2	Vortex 20/21 Ga
Description	RIS 2.2 mm OD 0.93 mm ID 0.72-0.50 mm	RIS 2.2 mm OD 0.94 (flute) / 0.84 mm (shaft) OD 0.71 (flute) / 0.58 mm (shaft)
30°	DP8730S	DP8930S
30° Angled		DP8930AS

Non-sterile. Limited re-use (10 cases)

MicroFlow® 2.2

30° **DP8730**

Required Sleeve: **DP5531**

24 per box **STERILE**
(12 pouches, 2 per pouch)

C-MICS Coaxial I/A Single Use

12 per box **STERILE**

85794ST

I/A handpiece curved shaft tip, with sleeve for 1.8 mm incision



85795ST

I/A handpiece 45° angled tip, with sleeve for 1.8 mm incision



2.8 mm Standard Incision Solution

Single Use Knives

E7602

1.8-2.2-2.8 mm
uni-trap



E7551A

2.65 mm,
angled up



E7559A

2.75 mm,
angled up



E7553A

2.85 mm,
angled up



E7515

15° sideport



6 per box **STERILE**

Capsulorhexis Forceps

8 per box **STERILE**

ET1202

Utrata Micro Capsulorhexis Forceps



ET2012

Devgan Micro Capsulorhexis Forceps




Phaco Pack Contents

6 per box **STERILE**

BL5110

Premium Vacuum Phaco Pack



1 VFM cassette assembly
and I/A set
1 Auxiliary drape
1 Screen drape
1 Remote control drape
2 Irrigation sleeves
1 Plug
1 Test chamber
1 Needle wrench

BL5120

Premium AFS Phaco Pack



1 AFM cassette assembly
and I/A set
1 Auxiliary drape
1 Screen drape
1 Remote control drape
2 Irrigation sleeves
1 Test chamber
1 Needle wrench

* the following phaco packs can also be used for 2.2 surgery:
BL5112, BL5111, BL5121

Phaco Needles

6 per box **STERILE**



Single use	19 Ga Standard	MicroFlow®	MicroFlow®Plus™	Vortex 19/20 Ga	Thin Tip 20 Ga*
Description	RIS 3.00-3.20 mm OD 1.70 mm ID 0.90 mm	RIS 2.65-2.75 mm OD 1.07 mm ID 0.90-0.51 mm	RIS 2.75-2.85 mm OD 1.17 mm ID 0.98-0.61 mm	RIS 2.65-2.85 mm OD 1.16 (flute) / 0.90 mm (shaft) ID 0.91 (flute) / 0.67 mm (shaft)	RIS 2.6-2.75 mm OD 0.90 mm ID 0.67 mm
30°	DP8130S	DP8230S	DP8330S	DP8830S	DP8430S
30° Angled				DP8830AS	DP8430AS
45°	DP8145S	DP8245S	DP8345S		DP8445S

Non-sterile. Limited re-use (10 cases)

19 Ga Standard	MicroFlow®	MicroFlow®Plus™	Thin Tip 20 Ga*
30° DP8130 45° DP8145	30° DP8230 30° Angled DP8230A 45° DP8245	30° DP8330 30° Angled DP8330A 45° DP8345	30° DP8430 30° Angled DP8430A 45° DP8445

* RIS-Recommended incision size in coaxial * OD-Outer diameter, * ID-Inner diameter

*Required Sleeve: **DP5511**

24 per box **STERILE** (12 pouches, 2 per pouch)

Coaxial I/A Single use

12 per box **STERILE**

85910ST

I/A Handpiece 45°
Silicone Coated Tip



85784ST

I/A handpiece 17 gauge
45° angled tip 0.3 mm
aspiration port



Biaxial I/A Single use

12 per box **STERILE**

85787ST

Irrigation handpiece 21 gauge
Curved shaft, smooth tip
2 lateral irrigation ports of 0.6 mm



85780S

Aspiration handpiece
21 gauge Curved shaft,
roughened tip



Stellaris accessories

Tray

BL3379
Stellaris® ultrasonic handpiece tray



Anterior Vitrectomy Packs

DP4400-6
20 gauge pneumatic anterior cutter with irrigation



MVS3201-6
20 gauge pneumatic anterior cutter without irrigation



6 per box **STERILE**

Maintenance

CX7120
Universal maintenance kit



BL4361
Dust cover

Bipolar

CX9430
Stellaris® bipolar cord lemo connector



D8200
Straight bipolar forceps 0.5 mm tip, with lemo connector



E7928
Bipolar scleral pencil, 45° bevelled tip, with lemo connector



12 per box **STERILE**

System Options

BL2111 Remote Control



Modules

BL2120 AFM module
BL2120KIT AFM module + cover plate
BL2310 VFM module
BL2310KIT VFM module + cover plate
BL2330 Ultrasound module (phaco & bipolar)
BL2300 Multi Media Center (MMC) module

Foot control

BL4390 Foot control battery (spare)
BL4391 Wall charger without electric adapter
BL4393 Battery charger cradle
BL4394 Foot control backup cable

Stellaris accessories

I/A Series

BL3161
I/A Tip Curved Silicone 0.4 mm



BL3162
I/A Tip Straight Silicone 0.4 mm



BL3163
I/A Tip 45 Degree Silicone 0.4 mm



BL3164
I/A Tip Straight Metal 0.4 mm



BL3165
I/A Tip Curved Metal 0.4 mm



BL3166
I/A Tip 45 Degree Metal 0.4 mm



BL3167
I/A Tip Straight Metal 0.3 mm



BL3168
I/A Tip 45 Degree Metal 0.3 mm



DP9744
Capsule Polisher Tip
45 Degree Metal 0.3 mm



DP9745
Capsule Polisher Tip
45 Degree Silicone 0.3 mm



DP9742
Adapter for Titanium
Handpiece



DP9739
Soft Tip I/A Straight



DP9740
Soft Tip I/A Curved



DP9741
Soft Tip I/A 45 Degree



Diamond Knives

E0130 1.8/2.2
MICS Tru-Size Diamond Knife 1.8-2.2 mm
With Integrated Incision Calibration System



E0130 1.8 mm
Profile Angled Diamond Keratome 1.8 mm



E0130 2.2 mm
Profile Angled Diamond Keratome 2.2 mm



E0130 2.65
Profile Angled Diamond Keratome 2.65 mm



E0130 2.8
Profile Angled Diamond Keratome 2.8 mm



E0130 3.0
Profile Angled Diamond Keratome 3.0 mm



Visit our Surgical Instruments website

Visit www.StorzEye.com for a complete listing of MICS and standard cataract surgery instrumentation from Bausch & Lomb