


# Sterilisation Trays

Sterilisation Trays.....206



# Sterilisation Trays

E7406

## Plastic Tray, 150 x 250 mm

Sterilisable plastic tray with silicone matting for supporting and protecting instruments. 150 x 250 mm. For 10-20 instruments. Overall height 25 mm (incl. floor nobs).



E7407

## Diamond Knife Sterilising Tray

Diamond knife tray with silicone inserts to securely hold the diamond knife. For 2 knives. 6cm wide. 15.87 cm length. 2.54 cm depth.



E7409

## Plastic Tray, 80 x 160 mm

Sterilisable plastic tray with silicone matting for supporting and protecting instruments. 80 x 160 mm. For 4-7 instruments. Overall height 25 mm (incl. floor nobs).



E7411

## Double Layer Microsurgical Instrument Tray

Size 254 x 154 x 49 mm. Capacity 20-30 instruments. Double layer sterilisable plastic tray with silicone matting on both layers for supporting and protecting instruments.

