



Vitreotomy Lenses



Vitrectomy Lenses

99500

Kilp Lens

Plan-concave silicone for routine vitrectomy. With a silicone ring of 4 mm height which juts out 2 mm inwards to fit the corneal curvature and for firm suction. Lens Ø 11 mm. Ring Ø 14 mm.



99501

Kilp Replacement Silicone Ring for Lens

Replacement ring for use with the Kilp vitrectomy lens 99500.



99501Y1

Clemens Lens Holder

For use with the Kilp vitrectomy lens 99500. C-shaped holder with 15 mm diameter. 45° angled. Overall length 67 mm.

