


Visidrape® Drapes

Full Body Drapes

Visidrape® full body drapes feature wider widths and a three-sided channel bag for maximum fluid capacity. Full body drapes can be utilized in all ophthalmic procedures.

581087 Visidrape Full Body Incise Drape


Overall size: 140 cm x 244 cm (55 in x 96 in)
Incise area: 9 cm x 11.5 cm (3 1/2 in x 4 1/2 in)
Capacity: 500 ml



10 per box

581088 Visidrape Full Body Aperture Drape


Overall size: 140 cm x 244 cm (55 in x 96 in)
Aperture area: 5 cm x 7.5 cm (2 in x 3 in)
Capacity: 500 ml
Tenting band forms to lift drape off patient's mouth and nose for enhanced comfort.
Includes: 2 - Attached fluid collection bags (left and right eye).
1 - Adhesive wick
1 - Adhesive wipe
1 - Corneal light shield
1 - Tenting band



10 per box

581715 Visidrape Full Body Incise Drape with Plastic Head and Channel Bag

Overall size: 165 cm x 254 cm (65 in x 100 in)
Incise area: 9 cm x 11 cm (3.5 in x 4 in)
Channel Bag: 62 cm x 40 cm (24 in x 15.75 in)
Capacity: 1000 ml



10 per box

581722 Visidrape Non-Woven Full Body Aperture Drape

Overall size: 165 cm x 254 cm (65 in x 100 in)
Aperture area: 6 cm x 4 cm (2.25 in x 1.5 in)




10 per box

Half Body Drapes

Visidrape 3/4 and 1/2 body drapes include an attached fluid collection bag and wick to remove blood and fluid from the surgical area. Non-woven drapes are lightweight and breathable for optimum patient comfort. They are low linting and moisture resistant to keep patient dry.

581084 Visidrape Non-Woven Incise Drape with Channel Bag


Overall size: 140 cm x 127 cm (55 in x 50 in)
 Incise area: 9 cm x 11.5 cm (3 1/2 in x 4 1/2 in)
 Fluid collection bag with channels: 14.6 cm x 25.4 cm (5 3/4 in x 10 in)
 Capacity: 500 ml



10 per box

581086 Visidrape Non-Woven Aperture Drape with Channel Bag


Overall size: 140 cm x 127 cm (55 in x 50 in)
 Aperture area: 5 cm x 7.5 cm (2 in x 3 in)
 Fluid collection bag with channels: 14.6 cm x 25.4 cm (5 3/4 in x 10 in)
 Capacity: 500 ml



10 per box

581721 Visidrape Osher Half Body Aperture Drape

Overall size: 193 cm x 185 cm (76 in x 73 in)
 Aperture area: 9.9 cm x 8.6 cm (3.9 in x 3.4 in)
 Fluid Collection Bag: 26.7 cm x 22.9 cm (10.5 in x 9 in)
 Capacity: 1000 ml
 Includes: 1 - Fluid collection bag
 1 - Eye Wick




10 per box

581723 Visidrape Non-Woven Half Body Incise Drape* with Channel Bag

Overall size: 165 cm x 183 cm (65 in x 72 in)
 Incise area: 9 cm x 11 cm (3.5 in x 4 in)
 Channel Bag: 62 cm x 40 cm (24 in x 15.75 in)
 Capacity: 1000 ml

*Without wick.




10 per box

Mid-Size Drapes

Visidrape mid-size drapes include an attached fluid collection bag and wick to remove blood and fluid from the surgical area. Non-woven drapes are lightweight and breathable for optimum patient comfort. They are low linting and moisture resistant to keep patient dry. Embossed plastic drapes are made of translucent non-glare material that are impermeable to moisture yet light weight to keep patient dry and comfortable.

581071 Visidrape Embossed Plastic Incise Drape with Fluid Collection Bag

Overall size: 101.5 cm x 127 cm (40 in x 50 in)
 Incise area: 9 cm x 11.5 cm (3 1/2 in x 4 1/2 in)
 Fluid collection bag: 15.2 cm x 16.5 cm (6 in x 6 1/2 in)
 Capacity: 250 ml




10 per box

Visidrape® Drapes

581078 Visidrape Embossed Plastic Aperture Drape with Fluid Collection Bag


Overall size: 101.5 cm x 127 cm (40 in x 50 in)
Aperture area: 5 cm x 7.5 cm (2 in x 3 in)
Fluid collection bag: 15.2 cm x 16.5 cm (6 in x 6 1/2 in)
Capacity: 250 ml



10 per box

581072 Visidrape Non-Woven Incise Drape with Fluid Collection Bag


Overall size: 101.5 cm x 127 cm (40 in x 50 in)
Incise area: 9 cm x 11.5 cm (3 1/2 in x 4 1/2 in)
Fluid collection bag: 15.2 cm x 16.5 cm (6 in x 6 1/2 in)
Capacity: 250 ml



10 per box

581079 Visidrape Non-Woven Aperture Drape with Fluid Collection Bag

Overall size: 101.5 cm x 127 cm (40 in x 50 in)
Aperture area: 5 cm x 7.5 cm (2 in x 3 in)
Fluid collection bag: 15.2 cm x 16.5 cm (6 in x 6 1/2 in)
Capacity: 250 ml




10 per box

Mini Drapes

Ideal for refractive or office procedures.

581095 Visidrape Mini Plastic Incise Drape


Overall size: 41 cm x 41 cm (16 in x 16 in)
Incise area: 9 cm x 11.5 cm (3 1/2 in x 4 1/2 in)
Adhesive backing.



10 per box

581096 Visidrape Mini Plastic Aperture Drape


Overall size: 41 cm x 41 cm (16 in x 16 in)
Aperture area: 5 cm x 7.5 cm (2 in x 3 in)
Adhesive backing.



10 per box

581700 Visidrape Incise


Overall size: 10.2 x 20.4 cm (4 x 8 in)
Adhesive area: 11 cm x 16 cm (4 in x 6 1/4 in)
Small incise provides complete coverage of lashes, lids and brows, reducing the risk of infection. Converts aperture to incise drape.



10 per box

581710 Visidrape Dual Refractive Drape

Overall size: 40 cm x 26 cm (16 in x 10 1/4 in)
Aperture area: 2.5 cm x 3.75 cm (1 in x 1 1/2 in)
US Patent 6345621.



10 per box