



# Eye Shields



# Eye Shields

99272

## Klein Eye Shield

Child size, white, 13 holes of 4 mm Ø. Non-sterile. Autoclavable, 45 x 55,5 mm.



99273

## Klein Eye Shield

For newborn, white, 13 holes of 3 mm Ø. Non-sterile. Autoclavable, 39,5 x 44,5 mm.



99275

## Klein Eye Shield

White, 13 holes of 6 mm Ø. Non-sterile. Autoclavable, 75 mm. Packaged 50/box.



99281Y4

## Klein Eye Shield

Transparent. With 20 mm Ø central hole. Non-sterile. 75 x 65 mm.



99441Y1

## Cornea Silicone Light Shield

For corneal protection, Non-sterile. Autoclavable. Blue, 18 mm.



99441Y2

## Cornea Silicone Light Shield

Transparent. For corneal protection. Non-sterile. Autoclavable. Light blue, 18 mm.



99441Y5

## Cornea Light Shield, black

For corneal protection. Non-sterile. Non thermo-stable.  
Do not heat above 60 °C. 18 mm.



E5699

## Crouch Corneal Protector Adult

Developed for oculoplastic procedures for corneal and scleral protection. Dimensions enable protector to rest on the sclera and not on the cornea. Easily removed with the E1500 Bishop-Harmon Forceps. Dimensions 23.5 mm x 25.8 mm. Individually packaged sterile. 10/box.



E5699 P

## Corneal Shields Pediatric

Developed by Dr. Earl Crouch for pediatric oculoplastic procedures for corneal and scleral protection. Dimensions enable protector to rest on the sclera and not the cornea. Easily removed with the E1500 Bishop-Harmon Forceps. Dimensions 17 mm diameter. Individually packaged sterile. 10/box.



ET4600

## Khan Scleral Shield

The Khan Shield manufactured in titanium helps to protect the globe during laser assisted blepharoplasty procedures. The oval shield was designed for patient comfort and includes anti-rotational posts and lash traction bar. Measures 23 mm x 28 mm.



