




Specula

Accessories	192
Locking/Adjustable	192
Wire	197



Specula

Accessories

93011

Spare Tubing For LASIK Specula

460 mm bifurcated silicone tubing. male luer-lock connector. Can be used for specula 52135, 52136, 52137, 52138, E4061, E4061 H, E4067.



93012

Aspiration specula tubing extension

For use with aspirating speculum Art. Nos. 52135, 52136, 52137, 52138. Silicone tube 2 m. With male and female luer connector.



99195

T-Connector Metal For Aspiration Speculum

For use with aspirating speculum Art. Nos. 52135, 52136, 52137, 52138.



E4118 T

Aspiration Tubing for Speculum

Replacement aspiration tubing for E4067. Packaged each.



Locking/Adjustable

52021

Weiss Eye Speculum

Child size. Right. Solid blades. 67 mm.



52022

Weiss Eye Speculum

Child size. Left. Solid blades. 67 mm.



Locking/Adjustable

52070

Sauer Pediatric Eye Speculum


Child size. Overall length 30 mm.



52111

Pretori Eye Speculum


Left. 68 mm.



52112

Pretori Eye Speculum

Right. 68 mm.



52131

Mellinger Eye Speculum

Overall length 70 mm.



52132

Mellinger Eye Speculum

Paediatric. 56 mm.



52135

LASIK Speculum with Aspiration

Angled 45°, open blades of closed wire loops. 17 mm blade outer width, 14.5 mm blade inner width, 33 mm max. Blade spread. 4 aspirating ports on lateral blade bends. Screw-type blade spread mechanism. 460 mm bifurcated silicone tubing. Male luer-lock connector.



Locking/Adjustable

52136

Speculum LASIK with aspiration, round valves

LASIK speculum with rounded valves for optimum patient comfort. Blade width allows for convenient insertion of suction ring. 3 aspirating ports on lateral blade bends. 17.5 mm blade outer length. 15 mm blade inner length. 38 mm maximum blade spread.



52137

Speculum LASIK with aspiration, open valves

LASIK speculum with rounded, open valves for optimum patient comfort. Blade width allows for convenient insertion of suction ring. 3 aspirating ports on lateral blade bends. 18 mm blade outer length. 15.5 mm blade inner length. 46 mm maximum blade spread.



52141

Mellinger/Beard Eye Speculum

Blade length 12 mm. Overall length 68 mm.



52142

Mellinger/Beard Eye Speculum

Child size. 53 mm.



52201

Cook Eye Speculum

9 mm blades. Small. 35 mm.



52202

Cook Eye Speculum

Medium. 38 mm.



Locking/Adjustable

52203

Cook Eye Speculum

Large, 48 mm.



52250

Bangerter Eye Speculum

55 mm.



SUH04

Adjustable Eye Speculum, sterile

Sterile single use adjustable eye speculum with 16mm solid blades. Packaged single. Sterile. 6/box.



52201Y3

Cook Eye Speculum

For infant, 5 mm blades. Small, 35 mm.



E4040

Castroviejo Speculum

Small, 15 mm blades, 28 mm blade spread, 75 mm.



E4041

Castroviejo Speculum

Medium, 16 mm blades, 31 mm blade spread, 82 mm.



Locking/Adjustable

E4050

Williams Eye Speculum 14 mm

Adult model. Spring-controlled spread of blades with locking mechanism. Dull finish. Blades 14 mm. Blade spread 52 mm. Overall length 86 mm.



E4051

Williams Eye Speculum 11 mm

Pediatric model. Spring-controlled spread of blades with locking mechanism. Dull finish. Blades 11 mm. Blade spread 41 mm. Overall length 69 mm.



E4056

Lancaster Eye Speculum

Spring-controlled spread of blades with locking mechanism. Dull finish. Blades 15 mm. Blade spread 41 mm. Overall length 69 mm.



E4064

Knapp Eye Speculum

Open 13 mm blades. Screw-controlled locking mechanism. Blade spread 60 mm. Dull finish. Overall length 85 mm.



E4065

Lieberman Speculum

Adjustable wire speculum with 17 mm open blades. Lid tension adjustable throughout procedure with thumb screw control. Straight shaft. Overall length 90 mm.



E4065 A

Lieberman Speculum Angled

Adjustable wire speculum with 17 mm open blades. Lid tension adjustable throughout procedure with thumb screw control. Angled 45°. Overall length 86 mm.



Locking/Adjustable

E4072

Mellinger Eye Speculum

16 mm blades. 27 mm blade spread. 70 mm.



x 1

E4073

Mellinger Eye Speculum

16 mm blades. 27 mm blade spread.
Self-locking mechanism. 70 mm.

x 1

E4094

Murdoch Eye Speculum

Fenestrated 14 mm blades. Self-locking mechanism.
Blade spread 24 mm. Dull finish. Overall length 57 mm.

x 1

E4111

Saunders Eye Speculum Premature Infant

Solid 4 mm blades. Locking thumb screw.
Blade spread 14 mm. Dull finish. Holes in bar decrease
weight of speculum. Overall length 35 mm.

x 1

Wire

SUH01

Barraquer Eye Speculum, sterile

Sterile single-use wire speculum with closed valves.
14mm blade length. Packaged single. Sterile. 24/box.


x24



52121

Barraquer Eye Speculum

14 mm blades. Overall length 40 mm.



x 1


Specula

Wire

52121H

Barraquer Eye Speculum Adult Heavy Wire

14 mm blades. 40 mm. As item 52121 but heavy valves
1.0 mm Ø.



52121K

Kratz-Barraquer Eye Speculum

14 mm blades. 40 mm. Open valves.



52123Y1

Barraquer-Osterhuis Eye Speculum


28 mm.



52125Y1

Barraquer Eye Speculum


Baby size. 5 mm blades. Blade distance 20 mm. 30 mm.



E4106KT

Temporal Approach Lid Speculum - Open Blade

Angled for temporal approach. Open 14 mm rounded
blades for patient comfort. Blade spread 18 mm.
Dull finish. Overall length 41 mm.



E4106N

Kratz-Nichamin Speculum

Open wire style 16 mm blade with a wing tip
twist giving the surgeon better on-axis accessibility.
Overall length 45 mm.



Wire

E4106 S

Barrquer Eye Speculum Solid Blade

Solid 14 mm blades. Blade spread 18 mm. Dull finish.
Overall length 41 mm.



E4106 ST

Temporal Approach Lid Speculum Solid Blade

Angled for temporal approach. Solid 14 mm blades.
Blade spread 18 mm. Dull finish. Overall length 41 mm.



E4106 T

Temporal Approach Lid Speculum Wire Blade


Angled for temporal approach. Wire 14 mm rounded
blades for patient comfort. Blade spread 18 mm.
Dull finish. Overall length 38 mm.



E4107

Barrquer Eye Speculum Pediatric 11 mm


Wire 11 mm blades. Blade spread 19 mm. Dull finish.
Overall length 34 mm.



E4107 K

Kratz-Barrquer Eye Speculum Pediatric

Open 10 mm blades. Blade spread 19 mm. Dull finish.
Overall length 35 mm.



E4107 P

Barrquer Eye Speculum Pediatric 5 mm

Closed wire 5 mm blade with a wing tip
twist giving the surgeon better on-axis accessibility.
Overall length 30 mm.



Specula

Wire

E4107 PS

Barraquer Eye Speculum Pediatric Solid Blade

Solid 6 mm solid blade speculum for newborn or pediatric use. Overall length 30 mm.



E4110

Sauer Eye Speculum Infant

Solid 8 mm blades. Blade spread 20 mm. Dull finish. Overall length 35 mm.



E4112

Alfonso Eye Speculum Newborn

Open 5 mm blades. Blade spread 27 mm. Dull finish. Overall length 45 mm.

