



# Fixation Rings



# Fixation Rings

E2424

## Josephberg Hilton Stabilizing Plate



The serrated edge of this reusable Stabilizing Plate applies pressure to the globe to add firmness and stability while inserting 25 Gauge or 23 Gauge trocar/cannulas for vitrectomy. The slot is aligned to measure an approximate distance of 3.5 mm from the limbus which is the targeted pars plana zone. The center opening of the plate is large enough to insert the trocar/cannulas at an appropriate angle while keeping the eye steady and firm. This instrument makes it easier for doctors to transition to the smaller gauge vitrectomy procedures and is adaptable to right or left hand use. Overall length 81 mm.



E4033

## McNeil Goldman Scleral & Blepharostat Ring



Adult Size. Inner diameter 17 mm.  
Outer diameter 24 mm.  
Loops retract lids. Polished finish.



E403412

## Flieringa Fixation Ring 12 mm



Stainless steel. Polished finish.  
Available in 1 mm increments from 12 mm to 22 mm.



E403413

## Flieringa Fixation Ring 13 mm



Stainless steel. Polished finish.  
Available in 1 mm increments from 12 mm to 22 mm.



E403414

## Flieringa Fixation Ring 14 mm



Stainless steel. Polished finish.  
Available in 1 mm increments from 12 mm to 22 mm.



E403415

## Flieringa Fixation Ring 15 mm



Stainless steel. Polished finish.  
Available in 1 mm increments from 12 mm to 22 mm.



# Fixation Rings

E403416

## Flieringa Fixation Ring 16 mm

Stainless steel. Polished finish.  
Available in 1 mm increments from 12 mm to 22 mm.



E403417

## Flieringa Fixation Ring 17 mm

Stainless steel. Polished finish.  
Available in 1 mm increments from 12 mm to 22 mm.



E403418

## Flieringa Fixation Ring 18 mm

Stainless steel. Polished finish.  
Available in 1 mm increments from 12 mm to 22 mm.



E403419

## Flieringa Fixation Ring 19 mm

Stainless steel. Polished finish.  
Available in 1 mm increments from 12 mm to 22 mm.



E403420

## Flieringa Fixation Ring 20 mm

Stainless steel. Polished finish.  
Available in 1 mm increments from 12 mm to 22 mm.



E403421

## Flieringa Fixation Ring 21 mm

Stainless steel. Polished finish.  
Available in 1 mm increments from 12 mm to 22 mm.



E403422

## Flieringa Fixation Ring 22 mm

Stainless steel. Polished finish.  
Available in 1 mm increments from 12 mm to 22 mm.



# Fixation Rings

E9012

## Nichamin Fixation Degree Gauge



Open 14 mm fixation ring with fine teeth to prevent damage to conjunctiva. Gauge is marked every  $10^{\circ}$  to determine the extent of the arc to be incised. Two broad marks align with steep axis. Overall length 112 mm.



ET9019

## Thornton Fine Fixation Ring 16mm



Open ring design for globe fixation during various ocular procedures. 16 mm in diameter with blunt fixation points. Swivel handle provides increased flexibility. Titanium. Overall length 108 mm.



ET9020

## Thornton Fine Fixation Ring 13mm



Open ring design for globe fixation during various ocular procedures. 13 mm in diameter with blunt fixation points. Swivel handle provides increased flexibility. Titanium. Overall length 107 mm.



ET9044

## Devgan-Seibel Fixation Ring



The off-set design of this instrument provides a level perspective independent of hand position and eases placement on the globe. Marks on the top of the ring are spaced  $60^{\circ}$  on either side of the center of the ring opening, serving as a guide for placement of a paracentesis  $60^{\circ}$  from the temporal clear corneal incision in coaxial phaco, or for placement of two micro incisions  $120^{\circ}$  apart in bi-manual phaco. The closed side of the fixation ring features marks that are  $30^{\circ}$  apart to serve as a guide for LRI nomograms. Overall length 117 mm.

